PANACEA ELETTRONICA

Art. 2429/V - 2430/V - 2431/V

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Dry incubator

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**Heissluft
Inkubator**


Produced by C.B.M. S.r.l.

(CEE 93/42)

① - MANUALE D'USO	pag.2
② - USER MANUAL	page.6
③ - MANUEL D'EMPLOI	pag.9
④ - BETRIEBSANLEITUNG	S.12

Dichiarazione di conformità – Declaration of conformity
Déclaration de conformité– Konformitätserklärung

Avvertenze:

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DRY INCUBATOR- USER MANUAL

Remarks:

- ◆ These instructions have to be carefully kept near the equipment, in order to be consulted by all users.
- ◆ In case of lack of voltage or lighting of the alarm pilot lamp the equipment stops the incubation cycle.
- ◆ We recommend not to lean any object on the equipment.
- ◆ Temperature oscillations of few °C, during the cycle, have to be considered normal.
- ◆ Always place one of the perforated grids on the lower plane (1).

Description of the equipment :

Referring to the figure of page 2 the incubator components are described here below:

- (A) Master switch
- (B) "Ready machine" yellow pilot lamp

This lamp is on when the incubator is ready and it is waiting the start of a cycle

- (C) Alarm red lamp 
This pilot lamp turns on in order to signal some possible anomalies and it is matched to a BIP sound. To set aside the alarm push the START button (see Alarms Table).
- (D) START push button
The function of this push button is to start a cycle, if the pilot lamp (B) is turned on. The **START** push button works also as RESET key, if it is pressed during an incubation cycle the cycle is cancelled. Moreover, if an alarm has intervened, by pressing the **START** push button, the alarm is cancelled.
- (E) Switch push button 
The function of this push button is to alternate, on the display (H), the visualization of temperature and incubation time.
- (F) Temperature pilot lamp 
When this pilot lamp is turned on it shows that the value read and/or set on the display (H) is referred to the temperature.
- (G) time pilot lamp $\frac{1}{2}$
When this pilot lamp is turned on it shows that the value read and/or set on the display (H) is referred to the incubation duration, heating excluded.
If it flashes it means that INFINITE time has been set.
- (H) display
Cycle parameters are visualized on the display : temperature inside the chamber, incubation time, set temperature.
In order to reset the cycle data, this must be stopped.
 - The temperature can be set between 30°C and 70°C,
 - Time can be set between 1 minute and 600 minutes or infinite time.

To set the INFINITE time the push button **+** must be kept pressed up to 600 and then it is necessary to wait few seconds.

When the writing "INF" appear on the display wait 3 seconds for automatic storage.

When INFINITE time is selected the Lamp Time (G) flashes.

To remove INFINITE time press the button **—**.

ATTENTION: AT INFINITE TIME THE EQUIPMENT CAN BE SWITCHED OFF ONLY MANUALLY.

- **(I) (L) keys**

The function of these keys is to modify the temperature setting and the incubation time, increasing or decreasing the value. By pressing the **+** key once, with the incubator waiting (B pilot lamp turned on), the set temperature value or the set incubation time is visualized on the display.

From this moment, by pressing the **+** or **—** keys it is possible to modify the value.

After reaching the desired value, wait 3 seconds for automatic acquisition.

If one of these keys is pressed during the cycle, the set values of time or temperature can be read on the display but it is not possible to modify them.

- **(M) Heating pilot lamp**

This lamp indicates the switching on of the heating elements of the equipment.

User's instructions

1. Connect the plug into the 230V outlet. Switch on the master (A).
The display shows a short self test, including the Firmware Version.
The yellow pilot lamp (B) is turned on in order to show that the machine is ready for the cycle.
2. Open the two sliding chimneys on both sides in order to make an easier air circulation inside the incubator.
3. Place the material to be treated inside the incubator and close the door (do not open the door during the cycle).
4. Set the desired temperature and time of incubation,
 - **Set the temperature :**
with the temperature pilot lamp on, press the **+** or **—** keys, or keeping the push button pressed, the increase or the decrease of the value becomes quick. After reaching the desired value, wait 3 seconds for the storing of the set value.
 - **Set the incubation time:**
press the switch push button **6** the timer pilot lamp turns on $\frac{1}{2}$. Repeat the previous operations, in order to visualize and modify the time. Wait 3 seconds for the storing of the set value.
5. After the parameters setting, start the incubation cycle by pressing the "START" push button. The yellow pilot lamp (B) turns off.
6. After some minute the incubator arrives at the set temperature and keeps it for the set time. The temperature measured inside the incubation chamber is visualized on the display (with a $\pm 2^{\circ}\text{C}$ tolerance) or, by switching the **6** push button, the remaining time to the end of the cycle.
7. At the end of the set time, the cycle automatically stops and the display shows the writing END matched to few BIP sounds, informing that the incubator has finished the cycle.

IMPORTANT: The cycle is considered valid only if, at the end, the writing “End” appears.

8. At the end of the cycle close the sliding chimney in case some material is left inside.
9. If the equipment is not used, switch off the master (A) or the push button START.

ATTENTION: AT INFINITE TIME THE EQUIPMENT CAN BE SWITCHED OFF ONLY MANUALLY.

ALARM TABLE

Display/Defect	Cause	Remedy
AL 1 – Alarm heating element broken	1) Excessive load 2) Heating element broken	1) Reduce the load and repeat the cycle 2) Replace the faulty heating element
AL 2 – Alarm Low Temperature	Heating element broken	Replace the faulty heating element
AL 3 – Alarm Overheating	CPU card in Error	1) Open immediately the door 2) Replace the faulty CPU card
AL 4 - Probe broken	CPU card in Error	Replace the faulty CPU card
AL 5 – Black Out	The Black Out	Repeat the cycle

Cleaning

ATTENTION: Before cleaning or screwing out any part of the machine make sure that this one is disconnected from the feeding network

To make a good cleaning use only neutral detergents and rinse

The gaskets and all plastic parts can be cleaned with a cotton flock soaked in alcohol.

Do not use sponges or abrasive materials which can damage the surface.

To clean do not use detergents, disinfectants or other liquids containing chlorides as they are corrosive for stainless steel; in case of contact please rinse immediately.

Warranty

This instrument is designed and built according to the Community Directive CEE 93/42 “Medical equipments” and C.B.M. is responsible for safety, reliability and performances provided that possible modifications or repairs are carried out by authorized people, the corresponding electric plant complies with IEC regulations and the equipment is used in compliance with these instruction for use.

C.B.M. guarantees its incubator for one (1) year, in case of manufacturing defects or failures of components, though not due, according to C.B.M. evaluation, to the operators' wrongful use of the apparatus.

The manufacturing company is released by any responsibilities for damages to people or things due to failure to earth the plant.

The warranty does not cover parts subject to wear and repairs or replacements carried out by personnel not authorised by C.B.M.

Details, descriptions and pictures are merely indicative.

C.B.M. S.r.l. reserves the right to carry out ameliorative modifications to its own products without any previous notice.

Disposal Of Electric And Electronic Refuse

Waste from Electrical and Electronic Equipment - Directive 2002/96/EC (WEEE)

The symbol (shown above) that appears on this product means that it should not be disposed of together with normal household refuse but should be treated separately.

Electric and electronic equipment may contain substances that are harmful to the environment and to the health of humans and should therefore be disposed of at special waste disposal sites.

There are harsh legal penalties in force against persons found disposing of such equipment in an irresponsible manner.

When the equipment is no longer of any use, it may taken back to the outlet from which you intend to purchase a replacement.

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