NADAL® hCG Pregnancy Test (test cassettes)

REF 142000/ 142000N-10/ 152000
142002/ 152002

DE Gebrauchsanweisung	2	CZ Návod k použití	37
EN Instructions for use	7	FI Käyttöohje	41
FR Instructions d'utilisation	12	DK Brugervejledning	46
ES Instrucciones de uso	17	NL Gebruiksaanwijzing	51
IT Istruzioni per l'uso	22	Symbols	59
PL Sposób użycia	27	Our Teams	60
PT Instruções de Utilização	32		

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Ref. 152002, Urin/Serum Test (25 mIU/ml), Probenmaterial: Urin

NADAL® hCG Pregnancy Test	EIA		
		Positiv	Negativ
	Positiv	130	0
	Negativ	0	178
Total		130	178
			308

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Ref. 152002, Urin/Serum Test (25 mIU/ml), Probenmaterial: Serum

NADAL® hCG Pregnancy Test	EIA		
		Positiv	Negativ
	Positiv	169	0
	Negativ	0	250
Total		169	250
			419

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14. Referenzen

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136

3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of hCG, Human Reproduction Update 12(6): 769-84
6. <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08pm>
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

Rev. 0, 2016-12-22 OM/JB

This package insert contains information on the following NADAL® hCG Pregnancy Tests:

Ref.	Cut-off-Level	Packaging
142000	10 mIU/ml in urine (early test)	Test cassettes, individually sealed
142000N-10	10 mIU/ml in urine (early test)	Test cassettes, in minigrip transport bag
152000	10 mIU/ml in urine or serum (early test)	Test cassettes, individually sealed
142002	25 mIU/ml in urine	Test cassettes, individually sealed
152002	25 mIU/ml in urine or serum	Test cassettes, individually sealed

1. Intended Use

The NADAL® hCG Pregnancy Test (142002/152002) or the NADAL® hCG Pregnancy Early Test (142000/142000N-10/152000) is a rapid, visual immunoassay for the qualitative, presumptive detection of human chorionic gonadotropin (hCG) in human urine or serum specimens. The test is intended for use as an aid in early detection of pregnancy and designed for professional use only.

2. Introduction and Clinical Significance

Human chorionic gonadotropin (hCG) is a glycoprotein hormone produced by the developing placenta shortly after fertilisation and can already be detected in the early stage of pregnancy in urine and serum.

In normal pregnancies, hCG can be detected in both serum and urine as early as 7 days after conception. The hCG level continues to rise very rapidly, frequently exceeding 50 mIU/ml by the first missed menstrual period, and culminates in the 10th-12th weeks of pregnancy in the range of 30,000-290,000 mIU/ml.

The early appearance of hCG in urine and blood soon after conception and its subsequent rapid rise in concentration makes it an ideal marker for the early detection of pregnancy.

Elevated hCG levels, comparable to the early stages of pregnancy may be associated with a hydatidiform mole or choriocarcinoma. For this reason, these possibilities should be ruled out before a pregnancy is diagnosed.

The NADAL® hCG Pregnancy Test is a rapid, visual immunoassay for the detection of hCG at the concentration of 10 mIU/ml or 25 mIU/ml in urine (142000/142000N-10/142002) and/or serum (152000/152002).

The test uses anti-hCG antibodies to detect hCG in the urine (or serum). The specificity of this immune reaction rules out cross-reactivity with structurally-related hormones (hFSH, hLH and hTSH) at above the maximum anticipated physiological levels.

3. Test Principle

The NADAL® hCG Pregnancy Test detects human chorionic gonadotropin through visual interpretation of colour development on the internal test strip. Anti-hCG antibodies are immobilised in the test line region of the membrane. Anti-mouse antibodies are immobilised in the control line region. During testing, the specimen reacts with anti-hCG antibodies conjugated to coloured particles and pre-coated onto the

conjugate pad of the membrane. The mixture then migrates along the membrane by capillary action and interacts with reagents on the membrane. If there is sufficient hCG in the specimen, a coloured line will appear in the test line region of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result.

The appearance of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 30 NADAL® hCG Pregnancy test cassettes (incl. disposable pipettes)
- 1 package insert

5. Additional Materials Required

- Specimen collection container
- Timer

6. Storage & Stability

Test kits should be stored at 2-30°C until the expiry date printed on the packaging. The test cassette should remain in the sealed foil pouch until use. Do not freeze tests.

Care should be taken to protect the components of the test kit from contamination. Do not use test cassettes if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, specimen collection containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Only urine samples (and/or serum samples for 152000/152002) and no other liquids should be used for testing.
- Bring samples to room temperature (15-30°C) prior to testing.
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible

pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).

- Humidity and high temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

Urine samples:

- Collect a urine specimen in a clean and dry plastic or glass specimen collection container which is free from any preservatives. From the first day on of the missed period the hCG levels are normally so high that urine samples collected at any time of the day may be used. However, excessive fluid intake should be avoided (dilution effect).
- If the test is used at an earlier point in time, the first morning urine should be used, as it generally contains the highest hCG concentration. For tests with a cut-off level of 25 mIU/ml, testing can be carried out starting from 2 days before the expected period. For tests with a cut-off level of 10 mIU/ml, the earliest point in time for testing is 4 days before the expected period. However, the risk of false negative results increases the earlier the testing is carried out (see also point 12 'Limitations').
- Specimens and tests should be brought to room temperature for the correct execution of the test procedure. Urine samples exhibiting visible precipitates should be centrifuged, filtered or allowed to settle to obtain clear supernatant for testing. Only clear supernatant should be used for testing.
- If testing is not to be performed immediately, urine specimens should be stored at 2-8°C for up to 48 hours. For long-term storage, specimens should be kept below -20°C.
- Prior to testing, frozen specimens must be completely thawed and mixed well. Specimens should not be frozen and thawed repeatedly.

Serum sample (only for 152000 and 152002):

- Collect a blood specimen aseptically by venipuncture.
- To avoid haemolysis, separate serum from cruent as soon as possible. Use only clear, non-haemolysed specimens. Specimens containing particulate matter should be clarified by centrifugation prior to assay.
- Icteric, lipemic, haemolysed and contaminated sera may lead to incorrect results.
- Ideally, testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum specimens should be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C.
- Bring specimens to room temperature prior to testing. Prior to testing, frozen specimens must be completely thawed and mixed well. Specimens should not be frozen and thawed repeatedly.

9. Test Procedure

Bring tests and specimens to room temperature (15-30°C) before testing.

1. Remove the test cassette from the sealed foil pouch and place it on a clean, level surface. Label the test cassette with a patient or control identification. For best results the assay should be performed without delay.
2. Holding the pipette vertically, add 3 drops of specimen (approximately 120 µL) directly into the specimen well (S) of the test cassette. Start the timer. If a micropipette is utilised, a new pipette tip should be used for every sample.

or

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

3. Wait for the coloured line(s) to appear.

Reading time: For tests which are used only with **urine samples** (142000/ 142000N-10/142002), the reading time is 3-5 minutes. Tests which are used with **urine or serum samples** (152000/ 152002) should be read 5 minutes after adding the sample, though samples with high hCG concentrations may show positive results with these tests in as little as 3 minutes.

For all tests, the reading time should not exceed 5 minutes.

3-5 min

5 min

Note: With a longer reading time, the risk of false positive results increases, as increasingly lower hCG concentrations are detected with increasing incubation time. As a result, in rare cases, also basal concentrations of non-pregnant women, especially menopausal women, may be detected.

10. Result Interpretation

Positive:

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the test line region (T).

Negative:

Only one coloured line appears in the control line region (C). No coloured line appears in the test line region (T).

Invalid:

The control line fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit

immediately and contact your local distributor.

Note:

- The colour intensity in the test line region (T) may vary depending on the concentration of the analyte present in the specimen. Therefore, any shade of colour in the test line region should be considered positive. Note that this is a qualitative test only and it cannot determine the concentration of the analyte in the specimen.
- Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.
- In the case of a strong test line (T), the control line (C) may be very pale. The test result should nevertheless be considered valid.
- In the case of a negative test result (possibly with strongly diluted samples) where pregnancy is nonetheless suspected, the test should be repeated after 48-72 hours or, in addition, a quantitative test should be performed, since the hCG concentration may not be high enough. The first morning urine contains the highest hCG concentration.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® hCG Pregnancy Tests are for professional *in-vitro* diagnostic use and should only be used for the qualitative detection of human chorionic gonadotropin (hCG).
- Very diluted urine specimens, exhibiting low specific gravity, may not contain representative levels of hCG. If pregnancy is still suspected after a negative test result, a first morning urine sample should be obtained 48-72 hours later and tested.
- If tests are used earlier than before the first day of a missed period, there is a risk of false negative results due to individual fluctuations of hCG concentrations. According to the literature, if 10 mIU/ml tests are used 4 days before the expected period, only about 53-74% of pregnancies are initially detected when using urine samples. If pregnancy is suspected, testing should therefore be carried out again with morning urine at a later point in time. Alternatively, the hCG concentration can be determined using a more sensitive quantitative method.
- In healthy men and healthy, fertile, non-pregnant women, no hCG can be normally detected using NADAL® hCG Pregnancy Tests. Basal hCG values increase in the menopause, so that in rare cases, especially with 10 mIU/ml tests, false positive results may occur.
- Very low levels of hCG (<50 mIU/mL) are present in urine or serum shortly after the implantation of a fertilised ovum in the womb. However, since a significant number of first trimester pregnancies miscarry for natural reasons, a test result that is weakly positive should be confirmed using other clinical methods. If necessary, the test should be carried out again with the first morning urine at a later point in time. The result of the retesting should be clearly positive.
- If, after a weak positive result, the result of the retesting using the first morning urine at a later point in time is negative, a spontaneous miscarriage might be the cause. Approximately 22% of clinically-undetected pregnancies end with an unnoticed miscarriage and are perceived as a slightly delayed period.
- A number of conditions other than pregnancy, including trophoblastic diseases, such as molar pregnancy, choriocarcinoma, testicular tumors, prostate cancer and lung cancer as well as medication containing hCG cause elevated levels of hCG which can be detected as positive by the test.
- Genito-urinary infections can lead to a false-positive hCG test result. An increased number of leucocytes, a basic pH, an increased amount of nitrite, erythrocytes and proteins in urine among others are indicative of genito-urinary infections. In case of suspicion, please verify this by means of urinalysis strips or specific swab tests. Very suspicious urine samples should be verified by laboratory testing. In the case of a positive result, no pregnancy should be diagnosed until the aforementioned possibilities are excluded.
- The tests are calibrated for intact hCG. The proportions of intact hCG and other hCG forms in serum and urine may differ according to the pregnancy stage. Especially at the beginning of pregnancy, when hCG values are still low, this can lead, in rare cases, to deviating results for both sample materials. In rare cases, at the beginning of pregnancy, deviating results may also be observed with tests based on the detection of beta-hCG.
- For urine samples: In rare cases, at a more advanced stage of pregnancy (approximately week 5-8 of pregnancy), increased amounts of hCG-beta-core fragments may be present in the urine hindering the detection reaction and leading to false negative results in spite of high hCG concentration. [7]
- For serum samples: As with all assays that utilise anti-mouse antibodies, there is a possibility of interference with anti-mouse antibodies (HAMA) in the sample. Samples from patients who have received preparations of monoclonal antibodies for diagnosis or therapy may contain HAMA (immunologically interfering substances). Such samples may cause false positive or false negative test results. Heterophile antibodies in serum may also cause erroneous results.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but must only be made by the physician after all clinical and laboratory findings have been evaluated. It is therefore recommended to confirm the obtained results with another testing method or to interpret them in conjunction with clinical findings.

Expected values:

The reference range (serum, 97.5 percentile, beta hCG) of healthy, fertile, non-pregnant women is <3 mIU/ml. In postmenopausal women this value increases to <5 mIU/ml. However, some isolated higher values can be found in other literature. After fertilisation there is a sharp rise so that, approximately 10 days after fertilisation (4 days before the missed period), the values of ≥10 mIU/ml are often already reached. The doubling time varies during the course of the pregnancy and is approximately 2.5 days at the beginning of the pregnancy. Due to variations in liquid absorption there is a high individual variability for hCG concentrations in urine. In the course of the first trimester the concentrations of 30,000-290,000 are determined, which then decrease in the 2nd and 3rd trimester to values of 45,000-100,000 mIU/ml. After birth hCG serum values decrease and are no longer detectable after approximately 11-17 days. [1], [6]

13. Performance Characteristics**Detection limit and measuring range**

The lower detection limit of the NADAL® hCG Pregnancy tests is set at 10 mIU/ml (for 142000/142000N-10/152000) or 25 mIU/ml (for 142002/ 152002) (Calibration with the 4th IRP of hCG). NADAL® hCG Pregnancy tests show no prozone effect up to a maximum concentration of 600,000 mIU/ml.

Analytical specificity**Cross-reactivity**

In a cross-reactivity study with the structurally related hormones hLH (300 mIU/ml), hFSH (1000 mIU/ml) and hTSH (1 mIU/ml) the NADAL® hCG Pregnancy tests showed no cross-reactivity. The concentrations used are higher than the physiologically maximum expected values of healthy women.

Interference

hCG free urine samples and urine samples containing hCG concentrations at the detection limit were spiked with numerous compounds. None of the compounds interfered with result formation.

Substance	Concentration
Acetaminophen	0.2 mg/ml
Acetylsalicylic acid	0.2 mg/ml
Ascorbic acid	0.2 mg/ml
Atropine	0.2 mg/ml
Bilirubin	0.2 mg/ml
EDTA	0.2 mg/ml
Ethanol	1.0%
Caffeine	0.2 mg/ml
Gentisic acid	0.2 mg/ml
Glucose	20 mg/ml
Haemoglobin	10 µg/ml
Urea	20 mg/ml
Uric acid	0.2 mg/ml
Creatinine	0.2 mg/ml
Oxalic acid	0.4 mg/ml
THC	0.1 mg/ml

Diagnostic performance characteristics

In clinical studies NADAL® hCG Pregnancy tests were tested with clinical sample material in comparison with another rapid test with the identical cut-off or enzyme immunoassay (EIA). Relative diagnostic performance characteristics were

calculated on the basis of the results. The results are presented in the following tables.

Ref. 142002 (urine test 25 mIU/ml), sample: urine

		Other hCG rapid test		
		Positive	Negative	Total
NADAL® hCG Pregnancy Test	Positive	66	0	66
	Negative	0	142	142
	Total	66	142	208

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 142000/142000N-10 (urine test 10 mIU/ml), sample: urine

		Other hCG rapid test		
NADAL® hCG Pregnancy Test		Positive	Negative	Total
	Positive	79	0	79
	Negative	0	131	131
	Total	79	131	210

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 152002, urine/serum test (25 mIU/ml), sample: urine

		EIA		
NADAL® hCG Pregnancy Test		Positive	Negative	Total
	Positive	130	0	130
	Negative	0	178	178
	Total	130	178	308

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 152002, urine/serum test (25 mIU/ml), sample: serum

		EIA		
NADAL® hCG Pregnancy Test		Positive	Negative	Total
	Positive	169	0	169
	Negative	0	250	250
	Total	169	250	419

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 152000, urine/serum Test (10 mIU/ml), sample: urine

		EIA		
NADAL® hCG Pregnancy Test		Positive	Negative	Total
	Positive	94	0	94
	Negative	0	106	106
	Total	94	106	200

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

Ref. 152000, urine/serum test (10 mIU/ml), sample: serum

		EIA		
NADAL® hCG Pregnancy Test		Positive	Negative	Total
	Positive	94	0	94
	Negative	0	106	106
	Total	94	106	200

Relative diagnostic sensitivity: >99.9%

Relative diagnostic specificity: >99.9%

Overall agreement: >99.9%

14. References

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
6. <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08pm>
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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Concordance générale : >99,9%

14. Bibliographie

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589

5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of hCG, Human Reproduction Update 12(6): 769-84
6. <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08pm>
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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Concordancia general: >99,9%

14. Referencias

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
6. <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19-2:08-pm>
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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14. Bibliografia

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778

4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
6. <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08pm>
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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Celková shoda: >99.9%

14. Reference

- Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
- http://americanpregnancy.org/while-pregnant/hcg-levels/from 15-8-19, 2:08 pm
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
- Gronowski AM et al. (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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Samlet overensstemmelse: > 99,9 %

14. Referencer

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
6. <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08pm>
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCG β Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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Algehele overeenkomst: >99,9%

14. Referenties

- Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
- Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
- Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
- Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
- Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update 12(6): 769-84
- <http://americanpregnancy.org/while-pregnant/hcg-levels/from-15-8-19,2:08pm>
- Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
- Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

Rev. 0, 2016-12-22 JD

NADAL® hCG Pregnancy Test
(Ref. 142000/ 142000N-10/ 152000/ 142002/ 152002)

NADAL® hCG Pregnancy Test
(Ref. 142000/ 142000N-10/ 152000/ 142002/ 152002)

NADAL® hCG Pregnancy Test
(Ref. 142000/ 142000N-10/ 152000/ 142002/ 152002)

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consultense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico-diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegränsning	Temperaturbegrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchcode	Merkning
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il l ingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrekkelig for <n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel – UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH