AFINION™ LIPID PANEL

EN 5-14

FI 25-34

SV 45-54

DA 15-24

NO 35-44

REF 1116802

1116921 Rev. A 2019/01

SYMBOLS/ SYMBOLER/ SYMBOLIT/ SYMBOLER/ SYMBOLER

Conformity to the European directive 98/79/EC
on *in vitro* diagnostic medical devices
Overensstemmelse med EU-direktivet 98/79/EC
om *in vitro* diagnostisk udstyr til medicinsk brug
Yhdenmukainen Euroopan direktiivin 98/79/EC
kanssa koskien *in vitro* diagnostisia lääkinnällisiä laitteita
Overensstemmelse med EU-direktiv 98/79/EC
for *in vitro* diagnostisk medisinsk utstyr
Överensstämmer EU-direktiv 98/79/EG
om medicintekniska produkter för *in vitro*-diagnostik

In vitro diagnostic medical device
In vitro diagnostisk medicinsk udstyr
In vitro diagnostinen lääkinnällinen laite
Medisinsk utstyr til *in vitro* diagnostikk
Medicinteknisk produkt för *in vitro*-diagnostik

Catalogue number
Katalognummer
Luettelonumero
Katalognummer
Artikelnummer

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Warning
Advarsel
Varoitus
Advarsel
Varning

EN AFINION™ LIPID PANEL

For use with the Alere Afinion™ AS100 Analyzer/Afinion™ 2 Analyzer.
For professional near-patient testing and laboratory use.

PRODUCT DESCRIPTION

Intended use

The Afinion™ Lipid Panel is an *in vitro* diagnostic test for quantitative determination of total cholesterol (Chol), high-density lipoprotein (HDL) cholesterol and triglycerides (Trig) in whole blood, serum and plasma. Values for low-density lipoprotein (LDL) cholesterol, non-HDL cholesterol and Chol/HDL ratio are calculated by the Afinion Analyzer.

Cholesterol measurements are used in the diagnosis and treatment of disorders involving excess or low cholesterol in the blood and lipid and lipoprotein metabolism disorders.

Summary and explanation of the test

Elevated cholesterol is a major cause of coronary heart disease (CHD) and an important cardiovascular risk factor. Clinical trials show that lipid lowering therapy reduces the risk for CHD and LDL is identified as the primary target of therapy^{1,3,5,6}.

The National Cholesterol Education Program's (NCEP's) updated recommendations for cholesterol testing and management are presented in the third report of Adult Treatment Panel (ATP III). This report recommends that a fasting lipoprotein profile should be obtained once every 5 years for all adults aged 20 years or older. A lipid profile consists of total cholesterol, HDL cholesterol, triglycerides and LDL cholesterol¹.

Strong epidemiological evidence links low levels of serum HDL cholesterol to increased CHD morbidity and mortality and thus low HDL-cholesterol level is strongly and inversely associated with risk for CHD¹. Elevated serum triglycerides are also associated with increased risk for CHD. In addition, elevated triglycerides are commonly associated with other lipid and non-lipid risk factors¹.

ATP III identifies non-HDL cholesterol (total cholesterol minus HDL cholesterol) as a secondary target of therapy in persons with high triglycerides (< 2.26 mmol/L (> 200 mg/dL)). The goal for non-HDL cholesterol in persons with high serum triglycerides can be set at 0.78 mmol/L (30 mg/dL) higher than that for LDL cholesterol^{1,6}.

The overall guidelines on Cardiovascular Disease (CVD) prevention in clinical practice strongly recommend modulating the intensity of the preventive interventions according to the total cardiovascular risk⁶.

Principle of the assay

Afinion Lipid Panel is a fully automated assay for quantitative determination of Chol, HDL and Trig in whole blood, serum and plasma. LDL, non-HDL and Chol/HDL are calculated by the Afinion Analyzer.

The Afinion Lipid Panel Test Cartridge contains all reagents necessary for determination of Chol, HDL and Trig in whole blood, serum and plasma. The sample material is collected using the sampling device integrated in the test cartridge. The test cartridge is then placed in the Afinion Analyzer. The analyzer detects the type of sample material by inspection of the sampling device at the beginning of the assay. The sample is then diluted. If the sample material is whole blood, the hematocrit will be measured in order to correct for the volume of red blood cells in the sample volume. Whole blood is lysed and hemoglobin is measured photometrical. Hematocrit is proportional to the hemoglobin concentration. The diluted whole blood sample will be filtered through a composite filter to separate blood cells from the plasma fraction. The filtered fraction will be used for the HDL, Chol and Trig measurements.

Total Cholesterol

Total Cholesterol is measured by an enzymatic colorimetric method. Esterified and free cholesterol are enzymatically converted into cholest-4-en-3-one and hydrogen peroxide. The hydrogen peroxide is used by hydrogen peroxidase to couple a phenol and 4-aminoantipyrin to a red quinone-imine dye. The colour intensity is directly proportional to the concentration of free and esterified cholesterol in the sample.

Triglycerides

Triglycerides are measured by an enzymatic colorimetric method. Triglycerides are enzymatically converted into glycerol by lipoprotein lipase. Glycerol is then further catalyzed in 2 steps to di-hydroxy-acetone-phosphate and hydrogen peroxide. The hydrogen peroxide then reacts with 4-aminophenazone and 4-chlorophenol under the action of peroxidase to form a red dyestuff. The color intensity is directly proportional to the concentration of triglycerides.

HDL cholesterol

In a first reaction, anti-human apolipoprotein B (apoB) antibody (R1) binds to apoB present on all lipoproteins but HDL (i.e. non-HDL). The antibody protects non-HDL from being degraded by pegylated cholesterol metabolizing enzymes in the second reaction (R2). In the R2 reaction free and esterified cholesterol of HDL are converted into cholest-4-en-3-one and hydrogen peroxide. The hydrogen peroxide is used by peroxidase to couple 4-aminoantipyrin to F-DAOS and forms a blue colour complex. The colour intensity is directly proportional to the concentration of free and esterified HDL cholesterol.

LDL cholesterol

NCEP recommends calculating LDL by use of the Friedwald formula²:

$$\text{LDL (mmol/L)} = \text{Chol} - \text{HDL} - \text{Trig}/2.2$$

$$\text{LDL (mg/dL)} = \text{Chol} - \text{HDL} - \text{Trig}/5$$

This equation is not valid for samples that have Trig above 4.52 mmol/L (400 mg/dL), in non-fasting specimen, or in patients with type III hyperlipoproteinemia.

non-HDL cholesterol

The sum of VLDL (very low density lipoprotein) + LDL is called non-HDL cholesterol. It is calculated routinely as total cholesterol minus HDL:
 $\text{non-HDL} = \text{Chol} - \text{HDL}$

Non-HDL cholesterol reflects the concentration of cholesterol within all lipoprotein particles currently considered atherogenic^{1,9}.

Chol/HDL ratio

Many studies show that the total cholesterol/HDL cholesterol ratio is a powerful predictor of CHD risk. This ratio reflects two powerful components of risk. High total cholesterol is a marker for atherogenic lipoproteins, whereas low HDL cholesterol correlates with the multiple risk factors of the metabolic syndrome and probably imparts some independent risk¹.

$$\text{Chol/HDL} = \text{Total Cholesterol} / \text{HDL Cholesterol}$$

Kit contents (per 15 tests unit)

15 Test cartridges packaged separately in foil pouches

1 Package insert

Materials required but not provided with the kit


Alere Afinion AS100 Analyzer/Afinion 2 Analyzer

Afinion Lipid Panel Control

Standard blood collection equipment

Description of the test cartridge

The main components of the test cartridge are the sampling device and the reagent container. The test cartridge has a handle, a barcode label with lot-specific information and an area for sample ID. See figure and table below.

Component	Function/composition
1 Sampling device	15 µL plastic sampling device to be filled with sample material. For collection of patient sample or control.
a. Closed position	
b. Lifted position	
2 Capillary	15 µL
3 Reaction wells	
a. Capillary wiper	Device for capillary wipe-off.
Filter tube	Device for separation of blood cells from plasma.
b. Filter	Device for separation of blood cells from plasma.
c. Trig reagent	Enzymes, Adenosine-tri-phosphate (ATP) and components for the colour reaction.
d. Chol reagent	Enzymes, phenol and components for the colour reaction.
e. HDL-R2 reagent	Enzymes and components for the colour reaction.
f. Dilution liquid	Buffer.
g. HDL-R1 reagent	Enzymes, antibody and components for the colour reaction.
h. Lysis reagent	Buffer with detergent.
4 Handle	For correct finger grip.
5 Barcode label	Contains assay- and lot-specific information for the Analyzer.
6 Optical reading area	Area for transmission measurement.
7 ID area	Space for written sample identification.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Use gloves.
- Do not use test cartridges after the expiry date or if the test cartridges have not been stored in accordance with recommendations.
- Do not use the test cartridge if the foil pouch or the test cartridge has been damaged.
- In case of leakage, avoid contact with eyes and skin. Wash with plenty of water.
- Do not reuse any part of the test cartridge.
- The used test cartridges, sampling equipment, patient samples and controls are potentially infectious. The test cartridges should be disposed of immediately after use. Proper handling and disposal methods should be followed in accordance with local or national regulations.

Warning

Chol reagent (120 µL)	Contains polidocanol Causes serious eye irritation
Dilution liquid (230 µL) HDL-R1 reagent (250 µL)	Contains isothiazolin ketone May cause an allergic skin reaction

STORAGE

Refrigerated storage (2-8°C)

- The Afinion Lipid Panel Test Cartridges are stable until the expiry date only when stored refrigerated in sealed foil pouches.
- Do not freeze.

Room temperature storage (15-25°C)

- The Afinion Lipid Panel Test Cartridges can be stored in unopened foil pouches at room temperature for 14 days.
- Note the date of removal from the refrigerator and the new expiration date on the kit container.
- It is, however, recommended to store the test cartridges in the refrigerator and remove a limited number of test cartridges at a time.

Opened foil pouch

- The test cartridge must be used within 10 minutes after opening the foil pouch.
- Avoid exposure to direct sunlight.
- Avoid relative humidity above 90 %.

SAMPLE MATERIAL

Sample volume required for Afinion Lipid Panel is 15 µL.

The following sample materials can be used with the Afinion Lipid Panel Test:

- Capillary blood (from finger prick)
- Venous whole blood with anticoagulants (EDTA or Heparin)
- Serum
- Plasma (EDTA or Heparin)
- Afinion Lipid Panel Control

Important!

- When testing triglycerides, the person should fast for 9-12 hours before the sample is collected.
- If the testing opportunity is non-fasting, only the values for Chol, HDL, non-HDL and Chol/HDL will be usable.

Sample storage

- Capillary blood without anticoagulants cannot be stored.
- Venous whole blood with anticoagulants (EDTA or Heparin) may be stored up to 4 hours at room temperature (15-25°C) or refrigerated (2-8°C) for 4 days. Do not freeze.
- Serum and plasma can be stored refrigerated for 10 days.
- Serum and plasma can be stored frozen for 9 months if tubes are sealed properly.
- Consult the Afinion Lipid Panel Control Package Insert for storage of the control materials.

PREPARING FOR ANALYSIS

 Consult the Afinion Analyzer User Manual for detailed instructions on how to analyze a patient or control sample. The Afinion Lipid Panel Quick Guide also provides an illustrated step-by-step procedure.

- The Afinion Lipid Panel Test Cartridge must reach an operating temperature of 18-30°C before use. Upon removal from refrigerated storage, leave the test cartridge in unopened foil pouch for approximately 15 minutes.

- Open the foil pouch just before use.
- Do not touch the test cartridge optical reading area. Hold the test cartridge by the handle.
- Label the test cartridge with patient or control ID. Use the dedicated ID area.
- While holding the test cartridge by the handle, turn it upside-down once and return to normal position before use.

COLLECTING A SAMPLE

Important!






- The patient should sit quietly for about 5 minutes before collecting the sample.
- Always use gloves.
- Do not use a test cartridge that has been accidentally dropped on the floor or lab bench after sample collection.

Fingerstick

- A warm hand and good blood flow from the puncture site are essential to draw a good capillary sample.
- Select a skin puncture site on one of the center fingers of either hand.
- Clean the finger twice using alcohol.
- Dry thoroughly with a gauze pad before pricking the finger.
- Use a lancet to prick the finger at the selected site. Properly dispose the lancet.
- Squeeze the finger gently to obtain a drop of blood and wipe away this first drop since this drop is likely to be contaminated with tissue fluids.
- Squeeze the finger gently again while holding it in a horizontal or slightly downward position until a second large drop of blood forms. Do not milk the finger. The puncture should provide a free-flowing drop of blood. **Excessive squeezing of the finger may cause an erroneous result.**
- Fill the capillary according to the procedure described below.
- Wipe off any excess blood from the finger and apply direct pressure to the wound site until the bleeding stops with a clean gauze pad.

Filling the capillary

- Remove the sampling device from the test cartridge.
- Fill the capillary; hold the sampling device at a slightly upward angle (see picture), bring the tip of the capillary just beneath the surface of the patient sample (a, d) or control material (e). **Do not tilt the sampling device downwards (b).** Be sure that the capillary is completely filled, see arrow (c).
- Avoid air bubbles and excess sample on the outside of the capillary. Do not wipe off the capillary.
- Immediately replace the sampling device **carefully** into the test cartridge.
- Once the capillary is filled with sample, the analysis of the test cartridge must start within 1 minute.

Sampling from a tube

- Patient samples stored refrigerated can be used without equilibration to room temperature.
- Mix the sample material well by inverting the tube carefully 8-10 times before collecting a sample.
- Fill the capillary according to the procedure described above.

Sampling from the AFINION™

Lipid Panel Control vial

- The control material can be used without equilibration to room temperature.
- Mix the control material well by inverting the vial 8-10 times before collecting a sample.
- Fill the capillary according to the procedure described above.

ANALYZING A SAMPLE

- Analyze the test cartridge following the procedure described in the Afinion Analyzer User Manual.
- The analysis time is 7-8 minutes.

INTERPRETATION OF RESULTS

Interpret the Afinion Lipid Panel test results with careful consideration to the patient’s medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, analyze the Afinion Lipid Panel Controls and retest the sample using a new Afinion Lipid Panel Test Cartridge. If the result is still questionable, confirm the result with another method. Analyze control materials frequently to verify the performance of the Afinion Analyzer System.

Measuring range

Two different measuring units are in use for reporting Lipid Panel test results. The Afinion Analyzer displays the results in mmol/L or mg/dL:

	Chol	HDL	Trig
[mmol/L]	2.59 – 12.95	0.39 – 2.59	0.51 – 7.35
[mg/dL]	100 – 500	15 – 100	45 – 650

 **Important!** Patient’s Lipid Panel values should be reported in units consistent with national recommendations. Please contact your local supplier if the national recommendation is unknown. Consult the Afinion Analyzer User Manual for instructions on how to change the Lipid Panel measuring unit.

If the Chol, HDL and/or Trig concentrations is outside the measuring range, LDL, non-HDL and Chol/HDL will be reported as above or below a value. In some cases no value for LDL, non-HDL and Chol/HDL will be calculated. See the table below.

Symbol	Cause/Explanation
Chol < 100 mg/dL Chol < 2.59 mmol/L	Total cholesterol concentration is below the measuring range.
Chol > 500 mg/dL Chol > 12.95 mmol/L	Total cholesterol concentration is above the measuring range.
HDL < 15 mg/dL HDL < 0.39 mmol/L	HDL cholesterol concentration is below the measuring range.
HDL > 100 mg/dL HDL > 2.59 mmol/L	HDL cholesterol concentration is above the measuring range.
HDL – – –	HDL concentration cannot be measured. Trig > 7.35 mmol/L (650 mg/dL).
Trig < 45 mg/dL Trig < 0.51 mmol/L	Triglyceride concentration is below the measuring range.
Trig > 650 mg/dL Trig > 7.35 mmol/L	Triglyceride concentration is above the measuring range.
LDL < value LDL > value	The Chol, HDL and/or Trig concentration is outside the measuring range. See displayed values.
LDL – – –	The Trig concentration is above 4.52 mmol/L (400 mg/dL) or LDL cannot be calculated as both the Chol, HDL and/or Trig is outside the measuring range.
non-HDL < value non-HDL > value	Chol and/or HDL concentration is outside the measuring range. See displayed values.
non-HDL – – –	The non-HDL cannot be calculated as both the Chol and HDL are outside the measuring range.
Chol/HDL < value Chol/HDL > value	The Chol and/or HDL concentration are outside the measuring range. See displayed values.
Chol/HDL – – –	The Chol/HDL cannot be calculated as both the Chol and HDL are outside the measuring range.

Standardization

Chol and HDL is traceable to the National Reference System for Cholesterol (NRS/CHOL). Trig is traceable to a Centers for Disease Control and Prevention (CDC) reference method.

Reference range

NCEP’s recommendations for cholesterol testing and management presented in the ATP III report describe the following classification for cholesterol and triglyceride testing:

LDL Cholesterol	[mmol/L]	[mg/dL]
Optimal	< 2.59	< 100
Near optimal/ above optimal	2.59 - 3.34	100-129
Borderline high	3.35 - 4.12	130-159
High	4.13 - 4.91	160-189
Very high	≥ 4.92	≥ 190
Total Cholesterol		
Desirable	< 5.18	< 200
Borderline high	5.18 - 6.19	200-239
High	≥ 6.20	≥ 240
HDL Cholesterol		
Low	< 1.04	< 40
High	≥ 1.55	≥ 60
Serum Triglycerides		
Normal	< 1.70	< 150
Borderline high	1.70 - 2.25	150-199
High	2.26 - 5.64	200-499
Very high	≥ 5.65	≥ 500

Hematocrit

When the sample material is whole blood, the hematocrit (Hct) is measured in order to correct for the volume of red blood cells in the sample volume (Hct range 20-60 %).

Interference

The substances listed below were tested for interference with Chol, HDL and Trig. No significant interference (< 10 %) was observed up to the following concentrations:

• Acetaminophen	200 mg/L
• Acetylsalicylic acid	1000 mg/L
• Acetylcysteine	1590 mg/L
• Ampicillin	1000 mg/L
• Ascorbic acid	6 mg/dL
• Atorvastatin	600 µg/L
• Bilirubin	20 mg/dL
• Calcium dobesilate	0.7 mg/dL
• Cefoxitin	2500 mg/L
• Cyclosporine A	5 mg/L
• Cyclosporine C	5 mg/L
• Fluvastatin	2.97 mg/L
• Hemoglobin (hemolysis)	0.5 g/dL
• Heparin	3000 U/L
• Ibuprofen	500 mg/L
• Intralipid	10000 mg/L
• Levodopa	15 mg/L
• Lovastatin	216 µg/L
• Metformin	40 mg/L
• Methyldopa	1.4 mg/dL
• Metronidazole	200 mg/L
• Pravastatin	7.32 mg/L
• Rifampicin	64.3 mg/L
• Simvastatin	80.4 µg/L
• Theophylline	100 mg/L
• Tetracycline	50 mg/L
• Anticoagulants (EDTA and heparin) at concentrations normally used in blood collection tubes do not interfere.	

Important! It is possible that other substances and/or factors not listed above may interfere with the test and cause false results.

Limitations of the test

- Do not analyze diluted samples with the Afinion Lipid Panel.
- Do not analyze hemolyzed or coagulated samples.
- Calcium dobesilate interferes with Afinion Lipid Panel at therapeutic levels and results in too low results for Chol, HDL and Trig⁷.
- Methyldopa concentrations above 1.4 mg/dL interferes with Afinion Lipid Panel and results in too low Trig results. This is above toxic levels of Methyldopa and there is no interference at therapeutic levels⁸.
- Levodopa concentrations above 15 mg/L may give too low HDL and Trig results. This is above the drug concentration at therapeutic level⁷.
- Acetylcysteine concentrations above 1590 mg/L may give too low Trig results. This concentration is above therapeutic levels, including the use of acetylcysteine in antidote treatment of acetaminophen intoxications⁸.
- Methylamino antipyrine (4-MAP), an active metabolite of the drug Metamizole, does cause interference with Afinion Lipid Panel at therapeutic levels and resulting in too low HDL and Trig values¹¹.
- If the Hct value is outside the range 20-60 %, no Lipid Panel test results will be reported and an information code will be displayed (see “Troubleshooting”). In these cases serum or plasma samples are recommended for Lipid Panel analysis.
- Hand creams and soaps with glycerol may cause falsely high triglyceride results.
- The triglyceride test measures triglycerides and free glycerol. Free glycerol usually is less than 0.11 mmol/L (10 mg/dL)^{4,10}.

QUALITY CONTROL

Quality control testing should be done to confirm that your Afinion Analyzer System is working properly and providing reliable results. Only when controls are used routinely and the values are within acceptable ranges can accurate results be assured for patient samples.

It is recommended to keep a permanent record of all quality control results. The Afinion Analyzer automatically stores the control results in a separate log. Consult the Afinion Analyzer User Manual.

Choosing control material

Afinion Lipid Panel Control from Abbott is recommended for routine quality control testing. Consult the Afinion Lipid Panel Control Package Insert.

If controls from another supplier are used, the precision must be determined and acceptable ranges for the Afinion Analyzer System must be established.

Frequency of control testing

Controls should be analyzed:

- anytime an unexpected test result is obtained.
- with each shipment of Afinion Lipid Panel Test Kits.
- with each new lot of Afinion Lipid Panel Test Kits.
- when training new operators in correct use of Afinion Lipid Panel and the Afinion Analyzer.
- in compliance with national or local regulations.

Verifying the control results

The measured value should be within the acceptable limits stated for the control material. Consult the Afinion Lipid Panel Control Package Insert.

If the result obtained for the control is outside the acceptable limits, make sure that:

- patient samples are not analyzed until control results are within acceptable limits.
- the control vial has not passed its expiry date.
- the control vial has not been used for more than 8 weeks.
- the control vial and Afinion Lipid Panel Test Cartridges have been stored according to recommendations.
- there is no evidence of bacterial or fungal contamination of the control vial.

Correct any procedural error and retest the control material.

If no procedural errors are detected:

- Retest the control material using a new control vial.
- Examine the laboratory's quality control record to investigate the frequency of control failures.
- Ensure that there is no trend in out-of-range quality control results.
- Patient results must be declared invalid when controls do not perform as expected. Contact your local supplier for advice before analyzing patient samples.

TROUBLESHOOTING

To ensure that correct Lipid Panel results are reported, the Afinion Analyzer performs optical, electronic and mechanical controls of the capillary, the test cartridge and all individual processing steps during the course of each analysis. When problems are detected by the built-in failsafe mechanisms, the analyzer terminates the test and displays an information code.

The table below contains Afinion Lipid Panel specific information codes. Consult the Afinion Analyzer User Manual for information codes not listed in this table.

Code #	Cause
101	Hematocrit below 20 %
102	Hematocrit above 60 %

PERFORMANCE CHARACTERISTICS

Method comparison

A method comparison was performed between the Afinion Lipid Panel used with the Alere Afinion AS100 Analyzer and two automated laboratory analyzers. The study included 91-94 patient samples. A linear regression analysis was performed and the data are summarized in Table 1.

Table 1: Method comparison.
y: Afinion Lipid Panel (capillary blood from finger prick) vs.
x: automated laboratory analyzers (serum)

Method	Analyte	Number of samples	Regression line (mmol/L)	Correlation coefficient (r)
1	Chol	94	$y = 1.08x - 0.27$	0.95
	HDL	91	$y = 1.02x - 0.02$	0.98
	Trig	93	$y = 1.07x - 0.05$	0.98
2	Chol	94	$y = 0.99x + 0.10$	0.94
	HDL	91	$y = 0.99x + 0.04$	0.97
	Trig	93	$y = 1.10x - 0.01$	0.98

Precision

Within-device (total) precision of Afinion Lipid Panel used with the Alere Afinion AS100 Analyzer was determined according to the CLSI Guideline EP5-A2. Afinion Lipid Panel Control C I (1), Control C II (2) and one serum sample (3) were assayed for 20 days and one Li-heparin whole blood sample (4) was assayed for 5 days. Precision data are summarized in Table 2.

Table 2: Within-device (total) precision. N=number of days, CV=Coefficient of Variation

Total Cholesterol			
Sample	N	Mean (mmol/L)	CV (%)
1	20	4.65	2.7
2	20	5.72	2.5
3	20	2.96	2.3
4	5	3.85	4.4
HDL Cholesterol			
Sample	N	Mean (mmol/L)	CV (%)
1	20	1.84	3.0
2	20	1.28	3.1
3	20	0.63	3.5
4	5	1.54	4.2
Triglycerides			
Sample	N	Mean (mmol/L)	CV (%)
1	20	1.97	2.7
2	20	3.34	2.7
3	20	0.75	4.0
4	5	1.25	5.7

Performance testing with the AFINION™ 2 analyzer

The performance of Afinion Lipid Panel obtained with the Afinion 2 Analyzer has been demonstrated to be equivalent to the performance obtained with the Alere Afinion AS100 Analyzer.

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1. National Cholesterol Education Program (NCEP). Third report of the expert panel on detection, evaluation, and treatment of high blood cholesterol in adults (Adult Treatment Panel III). NIH Pub. No. 02-5215. National Heart, Lung, and Blood Institute; 2005. 284 p. Available from: <http://www.nhlbi.nih.gov/guidelines/cholesterol/atp3full.pdf>
2. Bachorik PS, Ross JW, National Cholesterol Education Program (NCEP) recommendations for measurement of low-density lipoprotein cholesterol: Executive summary. Clin Chem 1995;41:1414-20.
3. Burtis C. A, Ashwood E. R., Tietz Fundamentals of Clinical Chemistry, 5th ed.
4. Burtis CA, Ashwood ER, Bruns DE, editors. Tietz textbook of clinical chemistry and molecular diagnostics, 4th ed. St Louis (MO): Elsevier Inc.; c2006. 2412 p. ISBN-13:978-0-7216-0189-2. ISBN-10: 0-7216-0189-8.
5. Perk J, et al. European guidelines on cardiovascular disease prevention in clinical practice. European Heart Journal 2012;33: 1635-1701.
6. Reiner Z, Catapano A. L, et al. ESC/EAS Guidelines for the management of dyslipidaemias. European Heart Journal 2011 (32):1769-1818
7. Breuer J. Report on the symposium "Drug effects in clinical chemistry methods". Eur J Clin Chem Clin Biochem 1996; 34(4):385-6.
8. Clinical and Laboratory Standards Institute. Interference testing in clinical chemistry; Approved guideline-2nd ed. CLSI document EP7-A2. Wayne (PA); Clinical and Laboratory Standards Institute; c2005, 120p. ISBN: 1-56238-584-4.
9. Brunzell JD, Davidson M, Furberg CD, Goldberg BG, Howard BV, Stein JH, Witztum JL. Lipoprotein management in patients with cardiometabolic risk. Consensus statement from the American Diabetes Association and the American College of Cardiology Foundation. Diabetes care 2008;31(4):811-22.
10. Jessen RH, Dass CJ, Eckfeldt JH, Do enzymatic analysis of serum triglycerides need blanking of free glycerol? Clin Chem 1990 (36): 1372-1375
11. V. Vlahov et al.: Pharmacokinetics of metamizol metabolites in healthy subjects after a single oral dose of metamizol sodium. Eur J Clin Pharmacol (1990) 38;61-65

 **Abbott Diagnostics Technologies AS**
Kjelsåsveien 161
P.O. Box 6863 Rodeløkka
NO-0504 Oslo, Norway
www.abbott.com/poct

ISO 13485 certified company

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.
1116921 Rev. A 2019/01