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Wandlager

Dr. Mach GmbH & Co. KG, Flossmannstrasse 28, 85560 Ebersberg, GERMANY

Tel. +49 (0)8092 2093 0, Fax +49 (0)8092 2093 50

www.dr-mach.de, info@dr-mach.de

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Tel. +49 (0)8092 2093 0, Fax +49 (0)8092 2093 50

www.dr-mach.de, info@dr-mach.de




1 Pictograms and advice in this document

	This symbol indicates possible sources of danger. Please also note the safety instructions and the hazard specification in the associated installation and operating instructions from the manufacturer of the support arm system.
	This symbol indicates a possible danger from electrical current. Please also note the safety instructions and the hazard specification in the associated installation and operating instructions from the manufacturer of the support arm system.
	This symbol indicates important information and tips about usage.

2 Pictograms and advice on the equipment and on the packaging

	This symbol indicates that the instructions for use must be followed.
	Serial number of the product
	Part number of the product
	Address of the product's manufacturer or distributor
	Date of manufacture and country of manufacture
	CE conformity symbol
	This symbol indicates that this is a medical device
	Unique device identifier (ID) of the product
	Information about disposal of the device

	This symbol indicates the maximum permissible total weight of the product including the end device
	Temperature range for transport and storage
	Atmospheric humidity for transport and storage
	Air pressure for transport and storage

3 General safety instructions

	Pay attention to the instructions for use when handling the product.
	To avoid the risk of electric shock, this product may only be connected to a supply network with a protective earth conductor.
	A primary-side ON/OFF switch to isolate the system from the supply network must be provided on site. The switch must meet the requirements of IEC 61058-1 for nominal voltage peaks of 4kV.
	This product is not designed for operation in environments enriched with oxygen.
	The product may only be used for the intended purpose. Otherwise, the manufacturer will not be liable for personal injury or damage to property.
	Changes to the product are prohibited and will invalidate the manufacturer's certificate of conformity and all warranty claims.
	Use only the mains units (or transformers) approved or supplied by the manufacturer. Non-observance will void the conformity of the product and release the manufacturer from any claims under warranty.
	Installation, maintenance and repair work may only be carried out by the manufacturer or by specially trained staff.
	Maintenance must be carried out on the light at least every two years

	Additional equipment that is connected to medical electrical equipment, must conform to the relevant IEC- or ISO standards (e.g. IEC 60950 or IEC 62368 for data processing equipment). Moreover, all configurations must meet the requirements for medical electrical systems (see Section 16 of the latest version of IEC 60601-1). Anyone who connects additional equipment to medical electrical equipment, is configuring a medical system and is therefore responsible for ensuring that the system meets the requirements for medical electrical systems. In case of doubt, contact your local representative or our technical customer services.
	When positioning the products, there is a risk of injury (e.g. crushing) and collisions with other objects (inventory) or walls.
	Parts that fall off could injure the patient or lead to an infection of the wound area!
	Do not remove the rating plate or the warning labels!

4 Further information

All information on operation, maintenance and disposal of the respective device and all technical data can be found in the instructions for use and in the online manual for the respective product (for download link, see instructions for use)

All information on installation, setting, operating, maintaining and disposing of the extension arm and spring arms can be found in the documentation from the manufacturer of the support arm system.


5 Mounting

5.1 Mounting preparation

	De-energize the on-site power supply and secure it against being switched on again.
	Before starting work, check whether the intended fastening point is suitable for supporting the entire intended installation. In any case, have the load-bearing capacity confirmed by a structural engineer in the form of static proof.
	All building law requirements must be met before assembly!
	If the anchoring is improperly fastened, the stability of the entire installation can be impaired! There is a risk of fatal injuries due to the entire suspension falling if the fastening of the ceiling attachment fails!

The wall bracket should be installed at a height that ensures a free passage of **at least 200 cm** below the lowest point of the suspension.

The height of the permanent electrical connection (if present) is 70 mm below the upper hole, the distance between the holes is 190 mm.


The recommended mounting height of the wall bracket can be found in the following table:

Luminaire	Minimum height of the lower hole X	required room height * to guarantee the minimum passage height of 200 cm
Mach LED 120 (Swing suspension arm)	2274	2500
Mach LED 120 (AC2000/VALiA)	2110	2400
Mach LED 130 (Swing suspension arm)	2274	2500
Mach LED 130 (AC2000/VALiA)	2130	2520
Mach LED 150 (Swing suspension arm)	2080	2400
Mach LED 150 (AC2000/VALiA)	2080	2400
Mach LED 300 (until 2020)	2230	2850
Mach LED 2	2230	2850
Mach LED 3	2230	2850
Mach LED 6MC	2180	2800
Mach LED 300 (from 2021)	2180	2800
Instrumenten- / Geräteträger	2080	2700

* Uppermost point of the spring arm; Values for full cardanic suspension (if applicable)

5.2 Mounting of the Wall Bearing

Standard anchoring

- Mount the wall bearing with the dowels approved by the structural engineer (usually 12 mm in diameter). To do this, use the two marked fastening holes at a **distance of 190 mm**.
- Put on the cover **1** so that it engages in the locks.
- Fasten the cover **2** to the cover **1** with two screws **5**.
- On delivery, the cover **3** and the cover **2** are pre-assembled with the four fastening screws **4**.

6 Electrical connection

6.1 Preparing the electrical connection

All Dr. Mach luminaires with LED technology require an external power supply unit, which can either be mounted directly in the wall bracket or in an external housing (transformer box). The power supply unit is usually included in the scope of delivery for examination lights; it must be ordered separately for surgical lights.


In the case of the Mach LED 120, Mach LED 130 and Mach LED 150 lights, the power supply unit is usually installed in the wall bracket. The wall mount has a 2-meter-long connection cable with an earthing contact plug. The socket used must be installed in accordance with the requirements of IEC or VDE 0100-710.

With the lights Mach LED 2, Mach LED 3, Mach LED 6M and Mach LED 300, the power supply unit is usually supplied in a plastic housing (transformer box). This can be installed in an existing switch cabinet, connected to the on-site power supply and wired to the luminaire accordingly.


In the case of the LED 2SC, LED 2SC Hybrid, LED 3SC and LED 300 lights, the power supply unit can optionally also be installed in the wall bracket (this is not possible with all other OR light models). These lights are connected to the on-site power supply as with the examination lights (Mach LED 120, Mach LED 130 and Mach LED 150) via a 2-meter-long connection cable with a protective contact plug.

In order to guarantee the operating voltage specified in the instructions for use and in the manual for the lamp, the cross-section of the installation cable between the power supply unit and the lamp must be dimensioned accordingly. The maximum connection cross-section on the terminal strip in the wall mount is **1.5 mm²**. For this reason, the length of the cable between the power supply unit and the wall bracket **must not exceed 10 m**.

6.2 Connection of the lines

- Connect the power supply units in the wall storage (see illustration) on the primary side to the on-site power supply.
- Mount power packs in the transformer box, if necessary in an existing control cabinet, connect to the on-site power supply and then wire up to the wall mount of the luminaire (the maximum cable length is 10 m).

- To connect the PA line **1**, put the ring eyelet on the press-in threaded bolt **5** and secure it with the washer **2**, toothed washer **3** and nut **4** (follow the sequence shown)
- Connect the plug cable on the wall bracket to an existing isolated ground receptacle.
- If a permanent connection is required, the plug-in cable can be disconnected in the wall layer, the on-site power supply can be placed on the wall bracket through an empty conduit and connected to the primary side of the power supply unit.

Before the electrical connection of the transformer box or the power supply unit, a qualified electrician must check whether the building installation, including the cable diameter and fuses, has been carried out in accordance with national / international regulations.