SIEMENS

Chek-Stix® Positive and Negative Control Strips for Urinalysis

For use with Siemens Healthcare Diagnostics Reagent Strips and Tablets for Urinalysis

SUMMARY AND EXPLANATION / INTENDED USE: When used with Siemens Reagent Strips and Tablets for Urinalysis, Chek-Stix® Positive and Negative Control Strips provide defined results for glucose, bilirubin, ketone (acetoacetic acid), specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes.

Each Chek-Stix Positive Control Strip is a firm plastic strip to which are affixed seven separate analyte areas. These each contain one or more natural or synthetic ingredients that, when dissolved out of the analyte areas in a measured quantity of distilled or deionized^a water, produce a Positive Control Solution. Each Chek-Stix Negative Control Strip contains six separate analyte areas. The ingredients are dissolved out of the analyte areas to produce a Negative Control Solution.

Chek-Stix Control Solutions should be used on a regular basis as part of a urinalysis quality control program.

Major areas of use are:

- ✓ When used as a known positive or negative solution:
 - To determine if Siemens Reagent Strips and Tablets are reacting properly.
 - To confirm the user’s ability to properly perform and reliably interpret the reagent strip tests.
 - For demonstration and teaching.
- ✓ When used as a blind or unknown solution:
 - To develop proficiency in routine urinalysis performed in the laboratory from sample handling through test procedures to reporting of results.
 - To provide confidence in obtaining good results in routine urinalysis.

Read the insert carefully before use. This product is for professional, *in vitro* diagnostic use.

PROCEDURE:

Materials required, but not provided:

- Specimen tubes (URIN-TEK® or equivalent; approximately 16 x 100 mm)
- Distilled or deionized^a water

Directions:

1. Place 12 mL of distilled or deionized^a water in an appropriately labeled specimen tube. Do not use tap water.
2. Remove a Chek-Stix Control Strip from the bottle and replace the cap immediately and tightly.
3. Place the strip into the tube. Cap tightly.
4. Repeat Steps 1–3 if using a second control.
5. Gently invert the tube(s) back and forth for 2 minutes.
6. Allow the tube(s) to stand for 30 minutes at room temperature.
7. Invert one more time, then remove and discard the strip(s), according to your standard laboratory procedures.

NOTE: After reconstitution, the Control Solutions should be treated and tested in the same manner as a urine specimen. Refer to the direction insert for the Siemens Reagent Strips or Tablets being used for complete information. Do not dip more than 12 Siemens Reagent Strips into a single tube of control solution, as expected values may change.

RESULTS: Chek-Stix control solutions provide results with Siemens Reagent Strips and Tablets in the same manner as urine specimens. The expected results are given in the Table of Values.

IMPORTANT: Do not make any chemical additions to the control solutions. This could seriously affect results due to chemical interference and/or dilution.

STORAGE AND STABILITY: Store Chek-Stix Control Strips in the original, tightly capped bottle at temperatures between 15–30°C. Do not store the bottle in direct sunlight. Do not remove the desiccant from bottle. The Control Strips are stable until the expiration date shown on the bottle label.

Prepare the control solutions as instructed in the PROCEDURE section. The reconstituted solution is stable for 8 hours at 18–23°C (64–73°F), except for bilirubin in the Positive Control, which is stable for 3 hours. If your laboratory operates outside of this temperature range, Siemens recommends that the reconstituted solution be refrigerated at 2–8°C (34–46°F) to maintain the 8 hour stability. Allow refrigerated control solutions to equilibrate to ambient temperature prior to use. Test bilirubin within 3 hours of preparation (bilirubin has limited stability in aqueous solution; protection from light will help maintain its reactivity). Ketone reactivity in the Positive Control will increase with time, due to continued hydrolysis of the reactive ingredient.

LIMITATIONS OF PROCEDURE: Chek-Stix Control Strips are for use with Siemens Reagent Strips and Tablets for Urinalysis.

REACTIVE INGREDIENTS: (Exact percent composition varies from lot to lot dependent upon component reactivity)

Positive Control Strips: Glucose; crystalline bilirubin; sodium methylacetoacetate; bovine hemoglobin; bovine serum albumin; 5-(4-sulfobutoxy)-2-methylindole sodium salt; sodium nitrite; protease (fungal).

Negative Control Strips: Potassium phosphate monobasic and dibasic; sodium chloride; urea.

AVAILABILITY: Chek-Stix Control Strips are available in a Combo Pak ((REF)) 1364 (03943435), which contains 1 bottle each of Positive Control Strips and Negative Control Strips (25 strips/bottle). Chek-Stix Positive Control strips are also available separately as Product ((REF)) 1360 (01918778) (1 bottle of 25 strips). For further information, contact your local Chek-Stix distributor.

TEST	EXPECTED RESULTS WITH SIEMENS REAGENT STRIPS AND TABLETS	
	POSITIVE CONTROL	NEGATIVE CONTROL
Glucose GLU	100 mg/dL – 250 mg/dL (SI Units: 5.5 – 14 mmol/L)	Negative
Bilirubin BIL	≥ Small ^b	Negative
Ketone KET	≥ Trace	Negative
Specific Gravity SG	1.000 – 1.015 (adjusted for pH) ^c	1.010 – 1.025 (adjusted for pH) ^c
Blood BLO	Moderate – Large	Negative
pH pH	≥ 8.0 ^d	6.0 – 7.0
Protein PRO	Trace – 100 mg/dL (SI Units: Trace – 1.0 g/L)	Negative
Urobilinogen URO	≥ 2 mg/dL ^e (SI Units: ≥ 33 µmol/L) ^e	0.2 – 1.0 mg/dL (SI Units: 3.2 – 16 µmol/L)
Nitrite NIT	Positive	Negative
Leukocytes LEU	Trace – Moderate ^b	Negative
ICTOTEST®	Positive ^b	Negative

^a If using deionized water, it must meet the minimum quality specified in CLSI Guideline C03-A4 ‘Preparation and Testing of Reagent Water in the Clinical Laboratory.’

^b Color will be atypical due to the alkalinity, high concentration of bilirubin, and color of Chek-Stix control solution.

^c Value listed is for Siemens Reagent Strips **only**. Results determined using the T.S. meter are lower, due to the constituents in the Chek-Stix Control Strips and the manner in which the Control Strips are processed.

^d Actual pH of the reconstituted solution, when measured by pH meter, is approximately 10. Readings obtained with Siemens Reagent Strips will generally be the highest reading available for the visual or instrumental interpretation.

^e Expected results when read instrumentally. When read visually, the slightly atypical color may be interpreted to be 1 color block lower than is expected instrumentally.

If test results are questionable or inconsistent with expected findings, refer to the direction insert for the Siemens Reagent Strips or Tablets being used for complete information on testing procedures and troubleshooting.

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