154701 140726

HemoCue® Glucose 201 RT Microcuvettes

The HemoCue Glucose 201 RT Microcuvettes are designed for use with the HemoCue Glucose 201 RT Analyzer and the HemoCue Glucose 201 DM RT Analyzer (referred to as the HemoCue Glucose 201 RT Analyzer in this document) with or without plasma conversion. A HemoCue Glucose 201 RT Analyzer with plasma conversion multiply the measured whole blood glucose value by a factor of 1.11¹. Please read the relevant Operating Manual for proper use of the systems^{2,3}.

Intended Purpose/Intended Use

The HemoCue Glucose 201 RT system is used for quantitative determination of glucose in whole blood, supplementing the clinical evidence in the diagnosis and treatment of patients with diabetes as well as monitoring of neonatal blood glucose levels. The system should not be used on critically ill neonates in neonatal intensive care settings. The HemoCue Glucose 201 RT system is for *In Vitro* Diagnostic use only. The HemoCue Glucose 201 RT Analyzer is only to be used with HemoCue Glucose 201 RT Microcuvettes. For professional use only.

IVD Medical Device Directive

The HemoCue Glucose 201 RT system comply with the IVD Medical Device Directive 98/79/EC and carry the CE mark.

Principles of the method/procedure

Principle of the method

The chemical reaction in the microcuvette has two phases, hemolysis and the glucose reaction. The glucose reaction is a modified glucose dehydrogenase method in which a tetrazolium salt (MTT) is used to obtain a quantification of glucose in visible light. α -D-glucose is transformed to β -D-glucose using mutarotase. Glucose dehydrogenase acts as a catalyst for the oxidation of β -D-glucose, to form NADH, which in the presence of diaphorase produces a colored formazan with MTT.

Principle of the procedure

The system consists of an analyzer together with microcuvettes. The microcuvette serves both as a pipette and as a measuring cuvette and is for single-use only. A blood sample of approximately 4 μ L is drawn into the cavity by capillary action. The measurement takes place in the analyzer in which the transmittance is measured and the absorbance and glucose level calculated. The calibration of the HemoCue Glucose 201 RT system is traceable to the YSI 2300 STAT Plus method and ID GC-MS method⁴. The HemoCue Glucose 201 RT system is factory calibrated and need no further calibration.

Composition

The microcuvette is made of copolyester plastic and contains the following reagent: <75 μ g/g cuvette MTT (methylthiazolylidiphenyl tetrazolium), <130 μ g/g cuvette saponin, <40 μ g/g cuvette sodium fluoride, <525 μ g/g cuvette enzyme mix: mutarotase (bacterial), glucose dehydrogenase (bacterial), diaphorase (bacterial), NAD (nicotinamide-adenine-dinucleotide), and non-reactive ingredients.

Warning and precautions

In accordance with “Good Laboratory Practices”, intervention levels and normal values should be established by each laboratory when instituting a new method. The specimen collection should be done according to instruction for use to reduce the risk of incorrect results. In cases of severe hypotension and peripheral circulatory failure, glucose measurement from capillary samples may be misleading. In such circumstances it is recommended that the glucose level should be measured using venous or arterial whole blood⁵. The measurement needs to be started no later than 40 seconds after filling the microcuvette. If a second sample is to be taken, fill a new microcuvette from a new drop of blood. This should not be done until the measurement of the first sample is completed. Do not remeasure the microcuvette. If the results exceeds the measuring range of the system “HHH” will be displayed. For more information read the relevant Operating Manual^{2,3}. Always handle blood specimens with care as they may be infectious. Consult local environmental authorities for proper disposal

Storage and handling of the HemoCue Glucose 201 RT Microcuvettes

The microcuvettes are to be stored at 0–30 °C (32–86 °F). Unused HemoCue Glucose 201 RT Microcuvettes should be kept in the original package but once an individual package is opened the microcuvette must be used

immediately or discarded. Use the microcuvettes prior to the expiration date printed on the package.

Specimen collection and preparation

Capillary, venous or arterial blood may be used. To achieve optimal result quality the specimen collection technique is essential. Make sure, regardless of type of blood, to do collection according to instructions for use. If venous or arterial specimen are used one of the following anticoagulants EDTA, sodium heparin and lithium heparin with and without gel and the glycolysis inhibitors sodium fluoride, sodium oxalate and potassium oxalate shall be used. Glycolysis is a major concern in all glucose measurements. To minimize the effect of glycolysis, measure the blood sample as soon as possible after taking the sample. Samples of blood collected in containers with recommended anticoagulants should be analyzed within 30 minutes. Mix all specimen tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8–10 times by hand.

Note! Lithium heparin with gel can not be used after centrifugation.

Note! Plasma specimen can not be used.

Procedure

The operating temperature of the HemoCue Glucose 201 RT system is 15-27 °C (59-80 °F). Please read the relevant Operating Manual for proper use of the systems^{2,3}.

Quality control

The HemoCue Glucose 201 RT Analyzer has an internal quality control, the “selftest”. Every time the analyzer is turned on, it will automatically verify the performance of the optronic unit of the analyzer.

Follow local guidelines regarding quality control procedures. If a quality control test is to be performed, GlucoTrol-NG (available at different levels) or GlucoLin from Eurotrol B.V. is recommended by HemoCue. For further information regarding quality controls contact HemoCue AB.

Results

The displayed range is 0–27.8 mmol/L (0–500 mg/dL) for a whole blood system and 0–31 mmol/L (0–560 mg/dL) for a plasma equivalent system. The method used in the HemoCue Glucose 201 RT system has been demonstrated to be linear between 0–27.8 mmol/L (0–500 mg/dL) for a whole blood system and 0–31 mmol/L (0–560 mg/dL) for a plasma equivalent system. If the results exceeds the measuring range of the system “HHH” will be displayed. If the patient assessment does not align with results from the HemoCue Glucose 201 RT system the test results should be verified against a reference method.

Table 1 Glucose 201 RT accuracy results for glucose concentration <3.3 mmol/L (<60 mg/dL) versus YSI 2300 STAT Plus, capillary sampling (Results in whole blood)

Within ± 0.28 mmol/L (± 5 mg/dL)	Within ± 0.33 mmol/L (± 6 mg/dL)	Within ± 0.56 mmol/L (± 10 mg/dL)	Within ± 0.83 mmol/L (± 15 mg/dL)
5/5 (100%)	5/5 (100%)	5/5 (100%)	5/5 (100%)

Table 2 Glucose 201 RT accuracy results for glucose concentration ≥3.3 mmol/L (≥60 mg/dL) versus YSI 2300 STAT Plus, capillary sampling (Results in whole blood)

Within ± 5 %	Within ± 10 %	Within ± 15 %	Within ± 20 %
150/217 (69%)	205/217 (94%)	216/217 (>99%)	217/217 (100%)

Table 3 Glucose 201 RT precision results in SD and CV (%), capillary and venous samples (Results in whole blood)

Glucose	No of samples	% of total	SD _{duplicates} mmol/L (mg/dL)		CV _{duplicates} %	
			YSI 2300 STAT Plus	Glucose 201 RT	YSI 2300 STAT Plus	Glucose 201 RT
0-11.1 mmol/L (0-200 mg/dL)	313	78.3	0.057 (1.02)	0.138 (2.49)	0.9	2.3
11.2-22.2 mmol/L (201-400 mg/dL)	74	18.5	0.128 (2.31)	0.267 (4.81)	0.8	1.7
22.3-27.8 mmol/L (401-500 mg/dL)	13	3.3	0.308 (5.55)	0.345 (6.22)	1.3	1.4
Total:	400	100	0.093 (1.67)	0.179 (3.22)	1.1	2.1

Limitations of the method/procedure

- a) The measured blood glucose value of, grossly lipemic specimens (intralipids >500 mg/L), turbid specimens, specimens containing >7% methemoglobin and high concentrations of glucosamine (>1.44 mmol/L, >26 mg/dL) may give incorrect results and should be interpreted with caution.
- b) At normal oral administration of xylose the HemoCue Glucose 201 RT system does not have any interference (Xylose <1.44 mmol/L <22 mg/dL) . However the high dosage used during a xylose absorption test may give incorrect results and should be interpreted with caution.
- c) It is verified according to CLSI Document EP7-A2 that the HemoCue Glucose 201 RT system does not have any interference at therapeutic levels for oral administration of Ascorbic Acid. Intravenous administration of Ascorbic Acid that leads to blood concentrations >10 mg/dL of Ascorbic Acid may give incorrect results and should be interpreted with caution.
- d) The HemoCue Glucose 201 RT system has no interference in blood specimens with hematocrit levels from 20 % to 70 %.
- e) A HemoCue Glucose 201 RT Analyzer with plasma conversion multiply the measured whole blood glucose value by a factor of 1.11 and display a plasma equivalent result¹. Since this factor is based on the relationship between plasma and whole blood at normal hematocrit, care should be taken when evaluating results in situations where the hematocrit may be extreme.
- f) The measurement needs to be started no later than 40 seconds after filling the microcuvette.

Bibliography

1. Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12, No 4.1
2. HemoCue Glucose 201 RT Operating Manual
3. HemoCue Glucose 201 DM RT Instructions for Use
4. Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794-800
5. Atkin et al, Annals of Internal Medicine, 1991; 114:12, 1020-1024

Symbols used

	Caution		Batch code		Use by (Year Month Day)
	<i>In vitro</i> diagnostic medical device		Catalogue number		Date of opening
	Do not reuse		Temperature limitation		Consult instruction for use
	CE mark		This product is covered by one or more patents		

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Literatur

1. Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12, No 4.1
2. Bedienungsanleitung des HemoCue Glucose 201 RT
3. Bedienungsanleitung des HemoCue Glucose 201 DM RT
4. Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800
5. Atkin et al, Annals of Internal Medicine, 1991; 114:12, 1020–1024

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Bibliografie

1. Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12, No 4.1
2. Gebruiksaanwijzing HemoCue Glucose 201 RT
3. Gebruikershandleiding HemoCue Glucose 201 DM RT
4. Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800
5. Atkin et al, Annals of Internal Medicine, 1991; 114:12, 1020–1024

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Referenser

1. Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12, No 4.1
2. HemoCue Glucose 201 RT bruksanvisning
3. HemoCue Glucose 201 DM RT bruksanvisning
4. Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800
5. Atkin et al, Annals of Internal Medicine, 1991, 114;12, 1020–1024

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Referencer

1. Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12, No 4.1
2. Brugsanvisning for HemoCue Glucose 201 RT
3. Brugsanvisning for HemoCue Glucose 201 DM RT
4. Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800
5. Atkin et al, Annals of Internal Medicine, 1991, 114:12, 1020–1024

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Referanser

1. Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12, No 4.1
2. HemoCue Glucose 201 RT Operating Manual
3. HemoCue Glucose 201 DM RT Instructions for Use
4. Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800
5. Atkin et al, Annals of Internal Medicine, 1991; 114:12, 1020–1024

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- Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12, No 4.1
- HemoCue Glucose 201 RT Operating Manual
- HemoCue Glucose 201 DM RT Instructions for Use
- Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800
- Atkin et al, Annals of Internal Medicine, 1991; 114:12, 1020–1024

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Bibliografia

1. Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12, No 4.1
2. Manuale d'uso di HemoCue Glucose 201 RT
3. Istruzioni per l'uso di HemoCue Glucose 201 DM RT
4. Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800
5. Atkin et al, Annals of Internal Medicine, 1991; 114:12, 1020–1024

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