cobas HbA1c Test

Hæmoglobin A1c

REF 08038694 190

▽ 10

cobas®

SYSTEM cobas b 101

Dansk

Anvendelse

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Referencer

- Goldstein DE, Little RR, Wiedmeyer HM, et al. Glycated hemoglobin: methodologies and clinical applications. Clin Chem 1986;32:B64-B70.
- Goldstein DE, Little RR, Lorenz RA, et al. Tests of glycemia in diabetes. Diabetes Care 1995;18:896-909.
- Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. Diabetes Care 1994;17:938-939.
- Santiago JV. Lessons from the diabetes control and complications trial. Diabetes 1993;42:1549-1554.
- Flückiger R, Mortensen HB. Review: glycated haemoglobins. J Chromatogr 1988;429:279-292.
- Bunn HF, Gabbay KH, Gallop PM. The glycosylation of hemoglobin: relevance to diabetes mellitus. Science 1978;200:21-27.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med 1993;329:977-986.
- UK Prospective Diabetes Study (UKPDS) group. Intensive blood glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). Lancet 1998;352:837-853.
- The International Expert Committee. International Expert Committee Report on the Role of the A1C Assay in the Diagnosis of Diabetes. Diabetes Care 2009, 32:1327-1334.
- WHO. Use of glycated haemoglobin (HbA1c) in the diagnosis of diabetes mellitus: abbreviated report of a WHO consultation. World Health Organization, Geneva. WHO Guidelines Approved by the Guidelines Review Committee. WHO/NMH/CHP/CPM/11.1. 2011:1-25.
- ADA. Diagnosis and classification of diabetes mellitus. Diabetes Care. Jan 2010;33 Suppl 1:S62-69.
- Ryden L, Grant PJ, Anker SD, et al. ESC Guidelines on diabetes, pre-diabetes, and cardiovascular diseases developed in collaboration with the EASD: the Task Force on diabetes, pre-diabetes, and cardiovascular diseases of the European Society of Cardiology (ESC) and developed in collaboration with the European Association for the Study of Diabetes (EASD). Eur Heart J. Oct 2013;34(39):3035-3087.
- IDF. Global guideline for type 2 diabetes. International Diabetes Federation. Brussels. Accessed in Oct-2015 under <https://www.idf.org/sites/default/files/IDF%20T2DM%20Guideline.pdf>. 2012:1-123.
- Ehehalt S, Gauger N, Blumenstock G, et al. Hemoglobin A1c is a reliable criterion for diagnosing type 1 diabetes in childhood and adolescence. Pediatr Diabetes. Nov 2010;11(7):446-449.
- Chiang JL, Kirkman MS, Laffel LM, Peters AL. Type 1 diabetes through the life span: a position statement of the American Diabetes Association. Diabetes Care. Jul 2014;37(7):2034-2054.
- ADA. Standards of medical care in diabetes-2015: summary of revisions. Diabetes Care. Jan 2015;38 Suppl:S4.
- Tapp RJ, Zimmet PZ, Harper CA, et al. Diagnostic thresholds for diabetes: the association of retinopathy and albuminuria with glycaemia. Diabetes Res Clin Pract. Sep 2006;73(3):315-321.
- Xin Z, Yuan MX, Li HX, et al. Evaluation for fasting and 2-hour glucose and HbA1c for diagnosing diabetes based on prevalence of retinopathy in a Chinese population. PLoS One. 2012;7(7):e40610.
- Kobold U, Jeppsson JO, Duellfer T, et al. Candidate reference methods for hemoglobin A1c based on peptide mapping. Clin Chem 1997;43:1944-1951.
- Jeppsson JO, Kobold U, Finke A, et al. Approved IFCC reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002;40:78-89.
- Consensus Statement on the Worldwide Standardization of the Hemoglobin A1c Measurement. American Diabetes Association, European Association for the Study of Diabetes, International Federation of Clinical Chemistry and Laboratory Medicine and International Diabetes Federation Consensus Committee. Diabetes Care 2007;30:2399-2400.
- Nathan DM, Kuenen J, Borg R, et al. Translating the A1C assay into estimated average glucose values. Diabetes Care 2008;31:1473-1478.
- Miedema K. Influence of hemoglobin variants on the determination of glycated hemoglobin. Klin Lab 1993;39:1029-1032.
- Chang J, Hoke C, Ettinger B, et al. Evaluation and Interference Study of Hemoglobin A1c Measured by Turbidimetric Inhibition Immunoassay. Am J Clin Pathol 1998;109(3):274-278.
- Rohlfing C, Connolly S, England J, et al. Effect of elevated fetal hemoglobin on HbA1c measurements: four common assay methods compared to the IFCC reference method. Clin Chem 2006;52 Suppl 6:A108.
- Sacks DB, Bruns DE, Goldstein DE, Maclaren NK, McDonald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem. Mar 2002;48(3):436-472.
- Behan KJ, Storey NM, Lee HK. Reporting variant hemoglobins discovered during hemoglobin A1c analysis - common practices in clinical laboratories. Clin Chim Acta. Aug 2009;406(1-2):124-128.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.

cobas HbA1c Test

Hæmoglobin A1c

- 30 Siegel E, Kellner M. The clinical practice recommendations of the German Diabetes Association. Exp Clin Endocrinol Diabetes. Jul 2014;122(7):383.
- 31 Bechel W. Roche internal Study Report "Performance evaluation cobas b 101 & HbA1c disc". Study CIM RD001325. 14-Dec-2012, 2012.
- 32 Brown A, Reynolds LR, Bruemmer D. Intensive glycemic control and cardiovascular disease: an update. Nat Rev Cardiol. Jul 2010;7(7):369-375.
- 33 American Diabetes Association. Standards of Medical Care in Diabetes-2010. Diabetes Care, 2010; 33 Suppl. 1: S11-S61.
- 34 American Diabetes Association. Standards of Medical Care in Diabetes-2011[J]. Diabetes Care, 2011, 34 Suppl. 1: S11-S61.

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© 2016, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim