

PRODUCT DATASHEET

SAP by Curas™, Fluid Absorber, Pad

Brand	: CURAS®
Category	: Nursing
Range	: SAP by Curas™
Product REF	: 09602
Product Description	: SAP by Curas™, Fluid Absorber, Pad, 24x13cm/ 9.5x5.1in (M), 500gsm, WardPak™



PRODUCT DESCRIPTION

The SAP Pad is a fluid absorbent pad. It consists of 5 layers of material including tissue paper and superabsorbent polymer (SAP) that enable to absorb fluid in seconds.

INTENDED USE

The SAP pad allows fast absorption of fluids such as liquids and body fluids by crystallizing them into gel. It ensures safe disposal and plays an important role in infection prevention systems.

FEATURES & BENEFITS

- SAP powder is embedded within the pad evenly.
- It converts fluids into gel for easy handling and disposal; improves hygiene performance of nurses.
- It can help to reduce odour and nosocomial infections by reducing free air particles.
- It has a high absorption rate; one SAP pad is able to absorb ~450 ml of human waste/body fluids.
- The pre-cut SAP pad is packed in WardPak™ (aluminium bag) to prolong the shelf life of the product and to provide convenience for users.
- To provide rapid and consistent absorption of body fluids and prevent any spillage of liquid.

PRODUCT SPECIFICATIONS

Ref	Dimensions (W x L)	Average Weight/Pc	Packing
09602	23.3 cm	13 g ± 1 g	WardPak™

Composition

- The pad contains tissue paper, virgin wood pulp, and SAP.
- SAP is made of sodium polyacrylate.

Country of Origin

- China

Sterilization

- Non-sterile

Shelf Life

- 5 years
- 6 months after opening

Selected Standards & Guidelines

- ISTA 1a, MSDS No.2613040033

Recommended Storage Condition

- Store the product at room temperature and avoid exposure to extreme temperatures.
- Keep the product at a cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If regular waste, dispose the used product together with its content into a pulp macerator, or treat as ordinary (domestic) incontinence waste. If a macerator is not available, incinerate or dispose the product into a landfill.
- If the content is contaminated or contagious, treat it as clinical waste and dispose it in accordance with the clinical waste protocol.
- Do not dispose used products into drain, wash basins or sinks.

Product Symbols



PACKAGING INFO

REF	EAN-13 Barcode	WardPak™ (Inner bag)	Carton	Gross Weight
09602	5051400096020	25 pcs/bag	500 pcs/carton	~9 kg

Components/Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	PET, PE, Aluminium	--	(√)	(√)		√	
Outer Carton	Brown corrugated board	--	(√)	(√)		(√)	√

√ = Highly recommended (√) = Recommended

ADDITIONAL INFO

- HS Code: 3906909090
- Documentation available upon request. Kindly contact: info@curas.com