

## PRODUCT DATASHEET

### Disposable Lid for 3L General Purpose Bowl WardPak™

|               |                                     |
|---------------|-------------------------------------|
| Brand         | : CURAS®                            |
| Category      | : Nursing                           |
| Range         | : Hygiene Pulp™                     |
| Product Group | : Hygiene Pulp™ Disposable Supports |
| REF           | : 08132                             |



#### PRODUCT DESCRIPTION

Disposable lid for 3L general purpose bowl. Suitable for disposal in a macerator. Packed in WardPak™ for easy dispensing. Non-sterile.

#### INTENDED USE

Disposable lid is used as a cover for disposable pulp 3L general purpose bowl to contain odours and prevent accidental spills. Single use.

#### FEATURES & BENEFITS

- Disposable and single use can reduce risk of cross-contamination and Hospital-Acquired Infections.
- Contains odours and prevent accidental spills.
- Allows for discreet transportation of used general purpose bowl.
- Made from environmentally friendly recycled material.
- Macerator friendly for convenient disposal.
- Die-cut to best fit Disposable 3L General Purpose Bowl(REF 08032) for improved hygiene.
- Integrated with perforation line for easy handling of the lid (one-hand operation).
- WardPak™ (inner bag) protects product from dust and contamination and allows easy dispensing at ward level.

#### PRODUCT SPECIFICATIONS

| REF   | Diameter | Average Weight/Pc | Packing  |
|-------|----------|-------------------|----------|
| 08132 | 23.1 cm  | 12 g ± 1 g        | WardPak™ |

#### Composition

- Lid is made of cardboard.

#### Country of Origin

- China

#### Sterilization

- Non-Sterile.
- Can be sterilized by ETO and Gamma radiation.

**Shelf Life**

- 8 years

**Selected Standards & Guidelines**

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 9001:2008

**Recommended Storage Condition**

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

**Proposed Disposal Method**

- If regular waste, dispose of used product with lid into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

**Product Symbols**



**PACKAGING INFO**

| REF   | EAN-13 Barcode | WardPak™      | Carton          | Gross Weight |
|-------|----------------|---------------|-----------------|--------------|
| 08132 | 5051400081323  | 100 pcs / bag | 700pcs / carton | ~9 kg        |

| Components / Material Properties |                   |              | Recommendations on Waste Handling |          |          |         |       |
|----------------------------------|-------------------|--------------|-----------------------------------|----------|----------|---------|-------|
| Packaging                        | Material          | Contains PVC | Incinerate                        | Landfill | Macerate | Recycle | Reuse |
| WardPak™                         | PE                | --           | √                                 | (√)      |          | √       |       |
| Outer Carton                     | Corrugated Carton | --           | (√)                               | (√)      |          | √       |       |
| Transport Pallet                 | Wood              | --           |                                   |          |          | (√)     | √     |

√ = Highly recommended      (√) = Recommended

**ADDITIONAL INFO**

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: [info@curas.com](mailto:info@curas.com)