


PRODUCT DATASHEET

Disposable 2L General Purpose Bowl w/Lid WardPak™

| | |
|---------------|-------------------------|
| Brand | : CURAS [®] |
| Category | : Nursing |
| Range | : Hygiene Pulp™ |
| Product Group | : General Purpose Bowls |
| REF | : 08033 |



PRODUCT DESCRIPTION

Disposable general purpose bowl made of pulp material with water resistance properties. 2L capacity. Suitable for disposal in a macerator. Lids are included for emergency use. Packed in WardPak™ for easy dispensing. Non-sterile.

INTENDED USE

2L disposable general purpose bowl can be used for holding liquid contents, vomitus, oral cleaning or holding drained urine from urine drainage bag/system. It can also be used as commode pan liner on selected commode chairs. Single use only.

FEATURES & BENEFITS

- An environmental friendly product that is made of 100% recycled newsprint paper.
- Disposable and single use can reduce risk of cross-contamination and Hospital-Acquired Infections.
- Water resistant, for a minimum period of 4 hours (applicable to general purpose bowl only).
- Integrated graduation marking for measurement indication.
- Cut-out rim provide convenience and comfort for vomiting user.
- Rounded edges on the base prevent cracking and increase wet strength.
- Wide and reinforced rim for easy holding and handling when the content is full.
- Canfit easily on selected commode chairs.
- Flat base design for stable placement.
- Reinforced base maintain material strength during usage and transportation with contents.
- Designed to be stackable for easy storage and reduces space.
- Disposable lid for discretion, and helps to contain odours and prevent accidental spills.
- Inner bag, WardPak™ protects product from dust and contamination.

PRODUCT SPECIFICATIONS

| REF | Capacity | Dimensions (D x H) | Average Weight/Pc | Packing |
|-------|----------|---------------------|-------------------|----------|
| 08033 | 2L | 31.5 Ø x 9.4 cm (H) | 55g ± 2g | WardPak™ |

Composition

- General purpose bowl is made of 100% recycled newsprint paper with water-based fluid resistant compound.
- Lid is made of cardboard.

Country of Origin

- China

Sterilization

- Non-Sterile.
- Can be sterilized by EtO and Gamma radiation.

Shelf Life

- General purpose bowl and lid: 8 years

Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 9001:2008
- PAS29:1999 Fluid resistance test for a minimum period of 4 hours (water at 35°C±3°C).
- PAS29:1999 Maceration test with a hospital sluice room type macerator, normal cycle.

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If regular waste, dispose of used product into a pulp macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

Product Symbols

PACKAGING INFO

| REF | EAN-13 Barcode | WardPak™ | Carton | Gross Weight |
|-------|----------------|--------------|---|--------------|
| 08033 | 5051400080333 | 50 pcs / bag | Bowl: 200 pcs / carton Lid: 5 pcs / carton | ~12kg |

| Components / Material Properties | | | Recommendations on Waste Handling | | | | |
|----------------------------------|-------------------|--------------|-----------------------------------|----------|----------|---------|-------|
| Packaging | Material | Contains PVC | Incinerate | Landfill | Macerate | Recycle | Reuse |
| WardPak™ | PE | -- | √ | (√) | | √ | |
| Outer Carton | Corrugated Carton | -- | (√) | (√) | | √ | |
| Transport Pallet | Wood | -- | | | | (√) | √ |

√= Highly recommended (√) = Recommended

ADDITIONAL INFO

- Curas[®] Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com