

## PRODUCT DATASHEET

### Disposable 1L General Purpose Bowl w/Lid WardPak™

Brand	: CURAS®
Category	: Nursing
Range	: Hygiene Pulp™
Product Group	: General Purpose Bowls
REF	: 08031



### PRODUCT DESCRIPTION

Disposable general purpose bowl made of pulp material with water resistance properties. 1L capacity. Suitable for disposal in a macerator. Lids are included for emergency use. Packed in WardPak™ for easy dispensing. Non-sterile.

### INTENDED USE

Disposable general purpose bowl for holding liquid contents, vomitus and various other purposes, e.g., oral hygiene and draining urine from urine drainage bag/system. Single use only.

### FEATURES & BENEFITS

- An environmental friendly product that is made of 100% recycled newsprint paper.
- Disposable therefore reduces risk of cross-contamination and Hospital-Acquired Infections.
- Water resistant, for a minimum period of 4 hours (applicable to general purpose bowl only).
- Integrated graduation marking for indication of approximate volume.
- Reinforced base for extra stability on flat surface.
- Wider and reinforced rim for easy holding and handling when full.
- Designed to be stackable for easy storage and reduces space.
- Inner bag, WardPak™ protects product from dust and contamination.

### PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (L x W x H)	Average Weight/Pc	Packing
08031	1L	22.8 ø x 7.6 cm (H)	32g ± 2g	WardPak™

### Composition

- General purpose bowl is made of 100% recycled newsprint paper with water-based fluid resistant compound.
- Lid is made of cardboard.

### Country of Origin

- China

### Sterilization

- Non-Sterile.

- Can be sterilized by EtO and Gamma radiation.

**Shelf Life**

- General purpose bowl and lid: 8 years

**Selected Standards & Guidelines**

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 9001:2008
- PAS29:1999 Fluid resistance test for a minimum period of 4 hours (water at 35°C±3°C).
- PAS29:1999 Maceration test with a hospital sluice room type macerator, normal cycle.

**Recommended Storage Condition**

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

**Proposed Disposal Method**

- If regular waste, dispose of used product into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

**Product Symbols**



**PACKAGING INFO**

REF	EAN-13 Barcode	WardPak™	Carton	Gross Weight
08031	5051400080319	100 pcs / bag	Bowl: 400 pcs / carton Lid: 5 pcs / carton	~13.4 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	PE	--	√	(√)		√	
Outer Carton	Corrugated Paper	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly Recommended                      (√) = Recommended

**ADDITIONAL INFO**

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: [info@curas.com](mailto:info@curas.com)