

Suprasorb® G Gel Wound Dressing


New quality

Impressive results.*

Examples of use

Case study 1:

Marked tissue damage following a fall¹



Case study 2:

Necrotic extravasation injury in a neonate²



i

The Suprasorb G Gel Wound Dressing is extremely effective in the treatment of nociceptive*** wound pain.^{3,4}



Sources:

- 1) Sian Fumarola. 21 days in the life of Edna. Poster Presentation: Tissue Viability Society Conference, Peterborough, April 2008.
- 2) Trudie Young. The use of an ionic hydrogel in a neonate with a necrotic extravasation injury. Poster Presentation: EWMA Conference, Lisbon, May 2008.
- 3) Sylvie Hampton. Comparing local (controlled) and national (uncontrolled) evaluations of ActiformCool. Sylvie Hampton. Poster Presentation: Wounds UK conference, Harrogate, November 2004.
- 4) Sylvie Hampton, Andy Kerr, Lynn King. Reduction of pain in painful wounds: a multi-centred audit. Poster Presentation: EWMA Conference, Stuttgart, September 2005.
- 5) Schmitz Michael, de Lange Stefanie, Abel Martin. Performance assessment by a hydrogel wound dressing and comparison to the institutions standard hydrogel dressing. Annual Meeting of the European Wound Management Association, (EWMA), Bremen, 11.–13. May 2016.
- 6) Donation von Suprasorb G Gel-Kompresse, nach interner Prüfmethode, in vitro Daten (*donation of Suprasorb G Gel Wound Dressing, using an internal test method, in vitro data*).
- 7) A small study in healing rates and symptom control using a new sheet hydrogel dressing. S. Hampton. Journal Wound Care, Vol 13, No 4, July 2004; 297-300 (20 wounds).
- * with autolytic debridement
- ** by relieving pain
- *** pain is caused by damage to tissue and irritation of exposed nerve endings
- **** Phenoxyethanol used as a preservative

Please note the new ordering information

Suprasorb® G Gel Wound Dressing

New quality

sterile, individually sealed

Original quality		New quality		
Size (cm)	REF	Size (cm) NEW	REF NEW	Shipping unit (pcs.) NEW
5 x 7.5	20 470	5 x 6.5	33 630	5/400
10 x 10	20 471	10 x 10	33 631	5/150
20 x 20	20 472	20 x 20	33 632	3/60

Please note the new ordering information.

Please follow the instructions for use.



Lohmann & Rauscher GmbH & Co. KG · Westerwaldstraße 4 · D-56579 Rengsdorf

Lohmann & Rauscher GmbH · Johann-Schorsch-Gasse 4 · A-1141 Wien



Suprasorb® G Gel Wound Dressing

New quality

Easy to use, impressive results.*



www.Lohmann-Rauscher.com

9506444 1016 e

Suprasorb® G Gel Wound Dressing

New quality

High moisture donation for effective autolytic debridement.

L&R's proven Suprasorb G Gel Wound Dressing is now available in a new version.

Suprasorb G Gel Wound Dressing supplies moisture, relieves pain^{3,4} and is very comfortable for patients to wear. Doctors and nursing staff favour this high-quality Wound dressing because of its ease of application and user-friendly handling.⁵

The hydrogel gently loosens and traps sloughs, necrotic tissue and debris found at the wound bed level.^{1,2}

Feature

Benefit

gently removes necrotic tissue and absorbs it into the gel^{1,2}




effective autolytic debridement

Feature

Benefit

provides moisture while relieving wound pain^{3,4,6}



the moist wound environment promotes granulation and epithelisation

better quality of life for patients**


Feature

Benefit

easy to mould^{5,7}



easy and safe to apply

 Suprasorb G Gel Wound Dressing is very comfortable for patients to wear.


Suprasorb® G Gel Wound Dressing

New quality


Gentle, yet effective on wounds.^{1,2}

Mode of action

The hydrogel layer of the Suprasorb G Gel Wound Dressing loosens necrotic tissue and debris while being gentle on the wound.^{1,2}



Structure




Hydrogel**** { high water content: approx. 70% acrylic polymers: approx. 30%

Transparent carrier film made from polyethylene

White application aid made from polyethylene

Blue carrier grid made from polyethylene

 Extensive studies attest to the cooling action and soothing effect of the Suprasorb G Gel Wound Dressing.^{3,4}



Suprasorb® G Gel Wound Dressing

New quality


Impressively easy to use.

Application


- Remove the sterile dressing from the packaging.
- Peel away one of the two white application aids. **A**
- Place the Suprasorb G Gel Wound Dressing on the wound and peel away the other application aid. **B**
- Secure the Suprasorb G Gel Wound Dressing with a suitable dressing (with a film dressing, e.g. Suprasorb F, an elastic conforming bandage, e.g. Mollelast, or as part of treatment for underlying diseases with compression bandages, e.g. Rosidal).




Indications




Dry to slightly exuding wounds




Superficial wounds



Chronic wounds,
e.g. lower leg ulcers, decubitus



Acute wounds,
e.g. scalds and first- to second-degree burns

 The Suprasorb G Gel Wound Dressing can be applied beyond the wound edge and cut to size for wounds if required.