ETHICON

a Johnson Johnson company



SUTURE

Description: Monofilament For smooth passage through tissue.

Synthetic For minimal tissue reaction.

Non-absorbable Provides permanent tensile strength

retention in tissue.

Color: Blue or clear.

Range: 10-0 to 1 (USP). Supplied as needled sutures, armed

with, for example, VISI-BLACK* needles, MULTIPASS* needles or CC-needles. Also available with pledgets.

PROLENE* sutures are also available as HEMO-SEAL*, with suture diameter to needle diameter ratio close to 1. HEMO-SEAL* is particularly useful in reducing needle hole

leakage (for example in synthetic vascular grafts).

Indications: PROLENE* sutures are intended for use in general soft

tissue approximation and/or ligation, including use in

cardiovascular, ophthalmic and neurosurgical procedures. Typical areas of use include:

Cardiovascular anastomoses

•Subcuticular sutures

Full details in the Instruction for Use included in every package.

