

# + ALLEVYN<sup>®</sup> ADHESIVE Foam Dressing

Trusted to stay in place  
with strong adhesion<sup>1-5</sup>



**Smith+Nephew**

ALLEVYN<sup>®</sup> ADHESIVE  
Foam Dressing

[www.smith-nephew.com/allevyn](http://www.smith-nephew.com/allevyn)

| EXUDATE LEVEL                                                                       |                                | PRIMARY INDICATIONS                                                                   |                                                                                       |                                                                                       |
|-------------------------------------------------------------------------------------|--------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |                                |  |  |  |
| MODERATE TO HIGH<br>ABSORBENCY                                                      |                                | LEG<br>ULCERS                                                                         | PRESSURE<br>ULCERS                                                                    | SURGICAL<br>WOUNDS                                                                    |
| DIABETIC<br>FOOT ULCERS                                                             | FIRST & SECOND<br>DEGREE BURNS | DONOR<br>SITES                                                                        | SKIN<br>TEARS                                                                         | FUNGATING<br>ULCERS                                                                   |
| OTHER INDICATIONS                                                                   |                                |                                                                                       |                                                                                       |                                                                                       |

Can be used in conjunction with **INTRASITE<sup>®</sup> GEL** for necrotic or sloughy wounds.

# ALLEVYN<sup>®</sup> ADHESIVE Dressings with our strongest adhesive and multi-way stretch are designed to stay in place.<sup>1-5</sup>

## Features and benefits

### ALLEVYN Range: absorbent, protective and gentle

ALLEVYN Foam Dressings are **absorbent, protective and gentle**, doing the basics you expect from a dressing and more.



#### Absorbent<sup>6-8</sup>

The hydrocellular foam core **innovatively balances absorption and retention** to create a moist wound healing environment<sup>3,9-12</sup> with minimised risk of maceration.<sup>\*3,9-11</sup>



#### Protective<sup>6,13</sup>

The breathable, waterproof<sup>6-8</sup> top film **provides a bacterial barrier<sup>14</sup>** while the dressing remains intact.



#### Gentle<sup>9,10</sup>

The wound contact adhesives are **easy to apply and remove<sup>3,11,15,16</sup>** and are shown to **minimise pain and trauma** during dressing changes and skin inspection.<sup>3,9</sup>

### ALLEVYN ADHESIVE



#### Secure adhesion

The low allergy/non-sensitising adhesive keeps the dressing **securely in place** without adhering to the wound.<sup>17-20</sup>



#### Multi-way stretch

Facilitates conformability and **helps to maximise patient comfort.**<sup>1,3,18</sup>



#### Cuttable<sup>\*21</sup>

Available in a wide range of sizes which can be cut to fit awkward areas such as elbows and knees.<sup>1,3,18,21,22</sup>

\*Cutting will compromise the bacterial barrier properties of the dressing.

†Cases, cartons and pouches compared to 2021 packaging dimensions.



**Less waste  
+ more care**



Our ALLEVYN ADHESIVE packaging is now 20% smaller.<sup>†23</sup>



## ALLEVYN<sup>®</sup> ADHESIVE Dressings



### ALLEVYN<sup>®</sup> ADHESIVE Foam Dressing

| Code     | Description              | Qty |
|----------|--------------------------|-----|
| 66000043 | 7.5cm x 7.5cm            | 10  |
| 66000348 | 7.5cm x 7.5cm            | 3   |
| 66000599 | 10cm x 10cm              | 10  |
| 66000044 | 12.5cm x 12.5cm          | 10  |
| 66000045 | 17.5cm x 17.5cm          | 10  |
| 66000744 | 12.5cm x 22.5cm          | 10  |
| 66000046 | 22.5cm x 22.5cm          | 10  |
| 66000700 | Small Sacrum 17cm x 17cm | 10  |
| 66000451 | Large Sacrum 22cm x 22cm | 10  |

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

**References:** 1. Amione P, Ricci E, Topo F. Comparison of Allevyn Adhesive and Biatain Adhesive in the management of pressure ulcers. *J Wound Care*. 2005;14(8):365–370. 2. Chamorro AM, Vidal Thomas MC, Mieras AS, Leiva A, Martínez MP, Hernández Yeste MMS. Multicenter randomized controlled trial comparing the effectiveness and safety of hydrocellular and hydrocolloid dressings for treatment of category II pressure ulcers in patients at primary and long-term care institutions. *Int J Nurs Stud*. 2019;94:179–185. 3. Leonard S, McCluskey P, Butters V, Winter R, Smith G. An evaluation of Allevyn™ Adhesive and Non-Adhesive foam dressings. *Wounds UK*. 2009;5(1):17–26. 4. Narci A, Embleton DB, Boyaci EO, Mingir S, Cetinkursun S. A practical offer for hypospadias dressing: Allevyn®. *Afr J Paediatr Surg*. 2011;8(3):272–274. 5. Peach V. Evaluating adhesive foam wound care dressings in clinical practice. *Wounds UK*. 2012;8(3):53–60. 6. S+N 2017. ALLEVYN Gentle Border Gen 2 - Physical Evaluation. Internal Report. DS/16/424/R V3. 7. S+N 2009. ALLEVYN Adhesive Dressing Physical Properties. Internal Report. DS/09/013/R7. 8. S+N 2007. Allevyn Ag Non-Adhesive Dressing Physical Properties. Internal Report. DOF 0702001. 9. Hurd T, Gregory L, Jones A, Brown S. A multi-centre in-market evaluation of ALLEVYN® Gentle Border. *Wounds UK*. 2009;5(3):32–44. 10. S+N 2022. Post Market Clinical Follow-Up (PMCF) Report for ALLEVYN GENTLE. CSD.AWM.22.050. 11. Rossington A, Drysdale K, Winter R. Clinical performance and positive impact on patient wellbeing of ALLEVYN Life. *Wounds UK*. 2013;9(4):91–95. 12. White R, Hartwell S, Brown S. Interim report on a study to assess the effectiveness and improved fluid uptake of new ALLEVYN. *Wounds UK*. 2007; 3(4):122–127. 13. Forri C, Gazineo D, Allegrini E, et al. Effectiveness of a multi-layer silicone-adhesive polyurethane foam dressing as prevention for sacral pressure ulcers in at-risk in-patients: Randomized controlled trial. *International Journal of Nursing Studies*. 2022;127:104172. 14. S+N 2023. Bacterial barrier properties of the ALLEVYN Dressing Range. Internal Report. CSD.AWM.23.032. 15. S+N 2017. An open, prospective, randomised, comparative volunteer trial to compare the performance of silicone adhesive dressings. Internal Report. GMCA-DOF-10. 16. S+N 2016. Product Performance of Next Generation ALLEVYN Life. Internal Report. (HVT080) GMCA-DOF/08. 17. Martini L, Reali U, Borgognoni L, Brandani P, Andriessen A. Comparison of two dressings in the management of partial-thickness donor sites. *J Wound Care*. 1999;8(9):457–460. 18. Bale S, Squires D, Varnon T, Walker A, Benbow M, Harding K. A comparison of two dressings in pressure sore management. *J Wound Care*. 1997;6(10):463–466. 19. Avanzi A, Martinelli M, Accardi S, Giraudi C, Peroli P. Adhesive Hydrocellular Dressing vs Hydrocolloid Dressing in the Management of 2<sup>nd</sup> and 3<sup>rd</sup> Degree Pus. Poster presented at: Eighth European Conference on Advances in Wound Management; April 26–28, 1998; Madrid, Spain. 20. S+N 2015. Product Safety Statement - ALLEVYN Adhesive. Internal Report. 21. S+N 2015. Cutting of ALLEVYN Variants. Internal Report. DS/14/318/R. 22. Williams C, Young T. Allevyn adhesive. *Br J Nurs*. 1996;5(11):691–693. 23. S+N 2022. Methodology for the reduction in packaging material of ALLEVYN® bordered dressings. Internal Report. CSD.AWM.22.045.