

+ PRIMAPORE[◇] Adhesive Non-Woven Wound Dressing is a low allergy adhesive¹ dressing for post-operative wounds

Designed to have a low adherent wound
contact layer and an absorbent pad²

Smith+Nephew

PRIMAPORE[◇]
Adhesive Non-Woven
Wound Dressing

<https://www.smith-nephew.com/en/>



Ask your pharmacist or doctor
about PRIMAPORE Dressing

PRIMAPORE[®] Adhesive Non-Woven Wound Dressing

is an absorbent post-operative dressing designed to help with the management of exudate²

Features and benefits



Designed to have a low adherent wound contact layer
PRIMAPORE Dressing is designed to have a low adherent wound contact layer for ease of removal



Low allergy adhesive¹
PRIMAPORE Dressing consists of a low allergy adhesive to help minimize skin reactions in patients^{1,3}



Conformable³
PRIMAPORE Dressing conforms to body contours.* The extensible non-woven layer is designed to help conformability of the dressing when in place⁴



Designed to help protect the wound
PRIMAPORE Dressing is designed to help protect the wound from external contamination



Breathable²
PRIMAPORE Dressing is breathable and permeable to moisture vapour²



Absorbent
PRIMAPORE Dressing is an absorbent dressing due to its pad²

Indications

PRIMAPORE Adhesive Non-Woven Wound Dressing is indicated to dress post-operative wounds; or cuts, lacerations or sutured wounds.



PRIMAPORE Adhesive Non-Woven Wound Dressing

S+N Code	Size	Carton
66000317	10cm x 8cm	20
66000318	15cm x 8cm	20
66000319	20cm x 10cm	20
66000320	25cm x 10cm	20
66000321	30cm x 10cm	20
66007140	35cm x 10cm	20
7133	7.2cm x 5cm	100
7135	8.3cm x 6cm	50

*As demonstrated *in vitro*.

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

Smith+Nephew
Kay Fiskers Plads 9
2300 København S
Denmark
T: +45 45 80 61 00

Smith+Nephew
Krokslätts Fabriker 39
431 37 Mölndal
Sweden
T: +46 31 74 65 800

Smith+Nephew
Snarøyveien 36
1364 Fornebu
Norway
T: +47 66 84 20 20

Smith+Nephew
Lentäjantie 1
01530 Vantaa
Finland
T: +358 (0)207 866 33

www.smith-nephew.com

®Trademark of Smith+Nephew
All Trademarks acknowledged
©2024 Smith+Nephew
43401-nord 06/24

References

1. Smith+Nephew 2018. Low allergy adhesive claims for CUTIPLAST Steril and CUTIPLAST Plus Steril. Internal report. PSS325. **2.** Smith+Nephew 2018. Non-woven Sample Manufacture Report for design verification testing and commercial evaluation. Internal report. C3188/R3. **3.** Smith+Nephew 2020. Claims Support - Low Risk of Skin Allergy and Irritation. Internal report. PSS 397. **4.** Smith+Nephew 2018. Verification testing of Usability related Design Inputs for Cutiplast (Sterile), Cutiplast Plus and Cicaplaie. Internal report. U/043/R2.