

+ IV3000[◇] Moisture Responsive Intravenous Catheter Dressings

Help prevent moisture accumulation, thereby minimising risk of infection and bacterial growth¹⁻³



Smith+Nephew

IV3000[◇]

Moisture Responsive Intravenous Catheter Dressing

www.smith-nephew.com

Intended use

- Frame Delivery
Suitable for dressing central line, peripheral, epidural and paediatric catheters.
- 1-Hand
Suitable for securing non-ported catheters, ported catheters and catheter stabilisation devices.

IV3000[®] Moisture Responsive Intravenous Catheter Dressings aid prevention of bacterial contamination,^{2,4-6} and have over 4 times higher moisture vapour transmission rate (MVTR) than other semi-permeable film dressings.^{7,8}

Features and benefits



Easy to use
IV3000[®] Dressings are easy to apply and remove whilst wearing surgical gloves^{9,10} and provide catheter security whilst being easy to remove.^{1,10,11}



Reduced dressing changes
A randomised, controlled trial demonstrated that use of IV3000 Dressings resulted in significantly fewer dressing changes.*¹²



Highly permeable
IV3000 Dressings shown to be significantly more permeable to water vapour than leading competitors.^{†13}



Minimised pain on removal^{10,14}
99% of patients reported no pain on dressing removal.^{‡10}



Reduced risk of infection
IV3000 Dressings maintain lower skin bacterial counts in comparison to other transparent film dressings.^{§2}

*n=25; p=0.03; compared to film dressing. †p<0.0001. ‡n=51.
§As demonstrated in a 12-volunteer study.

+ Evidence in focus

IV3000 Moisture Responsive Intravenous Catheter Dressing had the greatest MVTR and high conformability in a UK NHS clinical review of intravenous (IV) film dressings.

Key points



High conformability
with IV3000 Dressing (mean inflation pressure 130.5mmHg) versus most other dressings tested.⁷



6.5 times greater MVTR
with IV3000 Dressing than the mean for all dressings assessed (19,000 vs 2,928g/m²/24hr).⁷



>4.5 times greater MVTR with IV3000 Dressing versus the next best dressing (Tegaderm[™] IV. Advanced; 19,000 vs 4,102g/m²/24hr).⁷



Please, scan to see the full document

IV3000 Moisture Responsive Intravenous Catheter Dressings

S+N Code	Product	Size	Carton
4007	IV3000 1-Hand	6cm x 7cm	100
4008	IV3000 1-Hand	10cm x 12cm	50
4006	IV3000	7cm x 9cm	100
66004011	IV3000 Paediatric Dressing	5cm x 6cm	100
66004009	IV3000 Shaped	12cm x 9cm	50
4924	IV3000	6cm x 8.5cm	100
4973	IV3000	10cm x 14cm	50

S+N Code	Product	Size	Carton
4925	IV3000	10cm x 14cm	10
4649	IV3000	10cm x 20cm	50
4923	IV3000	6cm x 8cm	100
66800512	IV3000 Ported	11cm x 14cm	25
59410882	IV3000 Window Frame	10cm x 12cm	50
59410082	IV3000 Window Frame	6cm x 7cm	100

Smith+Nephew
Kay Fiskers Plads 9
2300 København S
Denmark
T: +45 45 80 61 00

Smith+Nephew
Krokslätts Fabriker 39
431 37 Mölndal
Sweden
T: +46 31 74 65 800

Smith+Nephew
Snarøyveien 36
1364 Fornebu
Norway
T: +47 66 84 20 20

Smith+Nephew
Lentäjantie 1
01530 Vantaa
Finland
T: +358 (0)207 866 333

®Trademark of Smith+Nephew
All Trademarks acknowledged
©2024 Smith+Nephew
42742-nord 04/23

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

References: 1. Keenlyside D. *Professional Nurse*. 1993;8(5):288–291. 2. Richardson MC. Paper presented at: Royal Society of Medicine Services International Congress and Symposium, Series No. 179; 1991; Harlow, Essex, UK. 3. Arrowsmith H. *British Journal of Nursing*. 1999;8(18):1231–1238. 4. Smith+Nephew 1990. Internal Report. SR/TW002/MB89-3. 5. Smith+Nephew 1994. IV3000 - Bacterial Barrier Claim. Internal Report. 6. Smith+Nephew 2005. Internal Report. DOF 0505005. 7. NHS Clinical Evaluation Team 2018. Clinical Review: Intravenous Vapour Permeable Film Dressings (IV Films) - Part One Securing Peripheral Cannulae in Adults (PVAD). UK. 8. NHS Clinical Evaluation Team 2018. Clinical Review: IV Film Dressings - Part Two Securing Central Venous Access Devices (CVAD). UK. 9. Smith+Nephew 2019. Internal Report. U/051/R1. 10. Wille JC, Blussé Van Oud Albas A, Thewessen EA. *The Journal of Hospital Infection*. 1993;23(2):113–121. 11. Maki DG, Stolz SS, Wheeler S, Mermel LA. *Critical Care and Medicine*. 1994;22(11):1729–1737. 12. Neufeld M. McMaster University: School of Graduate Studies, McMaster University; 1991. 13. Smith+Nephew 2023. Internal Report. CSD.AWM.23.005. 14. Mcie S, Pettitte T, Pridge L, Leeper D, Ostrow CL. *Am J Crit Care*. 2009;18(1):14–20.