



People.Health.Care.



Suprasorb H[®] Hydrocolloid-Dressing

Vedbæk, september 2019

Kære kunde,

Vi sender dette brev for at informere om, at vores Suprasorb H Hydrocolloid-Dressing får nye referencenumre og kommer i en forbedret produkt kvalitet. Den nye version er ydermere udvidet med tre indikationer. Produkt komposition, egenskaber og sterilisering er uændret.

Det vil være en løbende overgang fra den gamle version til den nye – startende fra oktober – I vil derfor kunne opleve at modtage tidligere referencenumre på nogle størrelserne.

Nye referencenumre

I bedes opdatere Jeres ordresystem med de nye referencenumre.

Produkt navn	Størrelse (cm)	Tidligere REF	Nye REF	Pakning (stk. i æske/karton)
Suprasorb H <i>Standard</i>	10x10	20403/121912	108830	10/660
Suprasorb H <i>Standard</i>	15x15	20401	108831	5/50
Suprasorb H <i>Standard</i>	20x20	20402	108832	5/50
Suprasorb H <i>Thin</i>	5x5	20424	108860	10/1650
Suprasorb H <i>Thin</i>	5x10	20410/121910	108861	10/480
Suprasorb H <i>Thin</i>	5x20	20411	108862	10/250
Suprasorb H <i>Thin</i>	10x10	20412/121911	108863	10/720
Suprasorb H <i>Thin</i>	15x15	20413	108864	5/50
Suprasorb H <i>Thin</i>	20x20	20414	108865	5/50
Suprasorb H <i>Border</i>	14x14	20420	108866	5/200
Suprasorb H <i>Sacrum</i>	14x16	20430	108867	5/200



People.Health.Care.



Indikationer Suprasorb H "NEW"

- ✓ Hudafskrabninger **NEW**
- ✓ Postoperative sår **NEW**
- ✓ Overfladiske forbrændinger **NEW**
- ✓ Tryksår
- ✓ Bensår
- ✓ Donorsteder
- ✓ Diabetiske sår

Forbedret produkt kvalitet

Suprasorb H "NEW" klarer sig bedre i test end Suprasorb H "OLD".

Test	Suprasorb® H * NEW	Suprasorb® H** OLD
Peel Adhesion 90° [N/cm]	3	2
FHC [g/m ² /24h]	3224	1897
Absorption [g/m ² /24h]	0,39	0,44
Thickness [mm]	5x5	20424

Source: R. Franke Prüfauftrag & Bericht; TR RD-D-2019-00130

*Suprasorb H new (thin 10x10)

** Suprasorb H old (thin 10x10)



People.Health.Care.

Suprasorb H "NEW" er uændret i produkt komposition, egenskaber og sterilisering

	Suprasorb H NEW	Suprasorb H OLD
Product composition	• Consist of a polyurethane backing coated with a hydrocolloid mass consisting of pectin, sodium carboxymethylcellulose, hydrocarbon resin, polyisobutylene, styrene-isoprene-styrene block copolymer, mineral oil and polybutene • The hydrocolloid layer is covered by an one-sided siliconized release paper	• Consist of a polyurethane backing coated with a hydrocolloid mass consisting of carboxymethylcellulose, polyisobutylene, paraffin, styrene isoprene-styrene-copolymer, resins, antioxidants and sorbitan ester • The hydrocolloid layer is covered by a protective paper siliconized on one side
Properties	• Self-adhesive wound dressing with semi-permeable carrier film	• Self-adhesive wound dressing with semi-permeable carrier film
Sterilisation	• Sterilized by radiation	• Sterilized by radiation

Spørgsmål

Hvis dette informationsbrev giver anledning til spørgsmål, er I meget velkomne til at kontakte os på mail: info@dk.LRmed.com eller på telefon: 45 89 02 25.

Med venlig hilsen
Lohmann & Rauscher A/S