



Wallace® Ring Pessaries

The Wallace® ring pessaries are available in a wide range of sizes to suit all patients. The flexible material provides comfort to the patient during fitting and duration of use.



Flexible PVC material

Sizes available (mm)

| | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|-----|-----|
| 50 | 53 | 56 | 59 | 62 | 65 | 68 | 71 | 74 | 77 | 80 | 85 | 90 | 95 | 100 | 110 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|-----|-----|

Designed to ease uterine prolapse and other uterine displacements.

Simple to Fit

- Ring softens when immersed in warm water to aid fitting in a figure of eight configuration

Easy to Maintain

- Remove and replace with a new ring pessary after three months

KEY PRODUCT BENEFITS

- Range of sizes to suit all patients
- Flexible PVC with large cross section
- Sixteen sizes from 50mm to 110mm (outside diameter)
- Individually packed in easy to open pouches

ORDERING INFORMATION

| Description | Quantity | Order Code |
|---|----------|--------------|
| Pessary Ring | 10 | 700/300/xxx* |
| Sizes available (mm) 50, 53, 56, 59, 62, 65, 68, 71, 74, 77, 80, 85, 90, 95, 100, 110 | | |
| *xxx is the size of the individual ring preceded by 0 if only two digits | | |

TECHNICAL INFORMATION

| | |
|-------------------|--------|
| Thickness of Ring | 12.9mm |
|-------------------|--------|