

## Storage & Transport HemoCue Microcuvettes

HemoCue microcuvettes shall be stored and transported at temperature conditions according to the table 1 below.

This is only valid for non-opened microcuvette packages at distributors and off-site storages and is not applicable for storage at end user customer sites.

Table 1. Storage and transport temperature conditions for HemoCue microcuvettes

| HemoCue microcuvettes     | Storage and Transport Temperature                                                                                                        | Microcuvettes can be stored or transported for a shorter period of time, maximum 2 weeks                                       |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Hb 201 and B-Hemoglobin   | +15 to +30 °C (+59 to +86 °F)                                                                                                            | -25 to +15 °C (-13 to +59 °F)                                                                                                  |
| Hb 301                    | +10 to +40 °C (+50 to +104 °F)                                                                                                           | -25 to +10 °C (-13 to +50 °F)                                                                                                  |
| Donor Hb Checker          | +15 to +30 °C (+59 to +86 °F)                                                                                                            | -25 to +15 °C (-13 to +59 °F)                                                                                                  |
| Plasma/Low Hb             | +15 to +30 °C (+59 to +86 °F)                                                                                                            | -25 to +15 °C (-13 to +59 °F)                                                                                                  |
| WBC and WBC DIFF          | +15 to +35 °C (+59 to +95 °F)                                                                                                            | -25 to +15 °C (-13 to +59 °F)                                                                                                  |
| Glucose 201 and B-Glucose | Below +8 °C (Below +46 °F)<br>It's recommended but not mandatory to store and distribute Glucose 201 and B-Glucose microcuvettes frozen. | Not applicable.<br>The microcuvettes shall always be stored according to labeling below +8 °C (+46 °F)                         |
| Glucose 201 RT            | +0 to +30 °C (+32 to +86 °F)                                                                                                             | -25 to +0 °C (-13 to +32 °F)                                                                                                   |
| Urine Albumin             | +2 to +8 °C (+35 to +46°F)                                                                                                               | Not applicable.<br>The microcuvettes shall always be stored according to labeling +2 to +8 °C (+35 to +46°F)<br>Do not freeze! |

If storage of microcuvettes is outside defined temperature limitations or specified time limit in table 1 above, customer service at HemoCue AB shall be contacted for further information.

Source: HCAB-000000640 [7]

2016.12.15