

neria™ soft

Intended use:

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency, apomorphine for Parkinson's Disease and morphine for pain management therapy. The neria™ soft infusion set has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion.


Product range:

Product code	Description	Cannula length	Tubing length	EAN	Pcs per inner box	Pcs per outer box
neria™ soft						
79-030-2522	neria™ soft infusion set	13 mm	30 cm	05705244000339	10	200
79-060-2522	neria™ soft infusion set	13 mm	60 cm	05705244000353	10	200
79-080-2522	neria™ soft infusion set	13 mm	80 cm	05705244000377	10	200
79-110-2522	neria™ soft infusion set	13 mm	110 cm	05705244000391	10	200
79-060-2622	neria™ soft infusion set	17 mm	60 cm	05705244000315	10	200
79-080-2622	neria™ soft infusion set	17 mm	80 cm	05705244001558	10	200
79-110-2622	neria™ soft infusion set	17 mm	110 cm	05705244001138	10	200
neria™ soft cannula						
76-000-2552	neria™ soft cannula	13 mm	N/A	05705244001091	10	200
76-000-2652	neria™ soft cannula	17 mm	N/A	05705244001534	10	200
neria™ soft tubing						
76-030-2652	neria™ soft tubing	N/A	30 cm	05705244001633	10	200
76-060-2652	neria™ soft tubing	N/A	60 cm	05705244001657	10	200
76-110-2652	neria™ soft tubing	N/A	110 cm	05705244001671	10	200

Product specifications:

Insertion angle 20-45 degrees
Soft cannula length 13 mm and 17 mm
Soft cannula outer diameter ~ 0,68 mm
Introducer needle 27 gauge
Tubing length in cm 30, 60, 80 and 110 cm
Tubing inner diameter 0,41 mm
Tubing outer diameter 1,50 mm

Material specifications:

PVC free Yes
Made with natural rubber latex No
Raw materials or derivation of animal origin No
The product contains phthalates No
Introducer & connector needles may contain . . ≤ 10% nickel

Certifications:

CE-marking Yes
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Tariff code HS9018390000
GMDN 35833

Sterilization:

Sterile Yes
Sterilized by Ethylene oxide
Shelf life 5 years
Usage time (drug depended) Maximum 72 hours

Drug compatibility/stability tested with:

Immunoglobulins
Apomorphine
Morphine (hydromorphone, morphine sulphate & morphine chloride)

Storage and handling:

	AVOID DIRECT SUNLIGHT - Keep the product away from direct sunlight and any humid environment
	KEEP DRY - Store in a dry place at room temperature
	DO NOT RE-USE
	DO NOT USE - If package is damaged
	CAUTION - See Instruction for Use



Unomedical a/s A ConvaTec Company

Aaholmvej 1-3
 Østved, 4320 Lejre, Denmark
www.neria.com