

PRODUCT DATA SHEET

INTRAFIX SAFESET LL,230CM(4063003)



Attribute Name	Attribute Value
Brand name	Intrafix®
Product name (HC)	Intrafix® SafeSet
Product type	Administration set
Application	Infusion
Type of administration set	Pressure
Container connector	IV standard spike
Length	230.00 cm
Inner-ø tube	3.00 mm
Outer-ø tube	4.10 mm
Type of filter	Fluid filter
Position of filter	Drip chamber
Pore size (Fluid filter)	15 µm
drops/ml	20
Air ventilation	vented

B. Braun Melsungen AG | 34212 Melsungen | Germany
Tel. +49 5661 71-0 | Fax +49 5561 71-4567 | www.bbraun.com

Attribute Name	Attribute Value
Flow control device	Roller clamp
with AirStop	Yes
with PrimeStop	Yes
Number of back check valves	-
Number of needle-free valves	-
Injection access	No
Connection	Luer-Lock/male
Light-protected	No
Color of light protection	-
Pump compatibility	Infusomat® P
max. pressure	2.0 bar
Material	PVC
Type of packaging	Peel pack
Product (unpack.) contains PVC	Yes
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product contains PVC	Yes
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Disposable / Reusable	disposable
Sterility	sterile
Type of Sterilization	E0