

Vitrex Flex 3

Depth Adjustable Safety Lancet

Technical Data Sheet



PRODUCT:

The Vitrex Flex 3 safety lancets are used to carry out skin punctures to collect micro-samples of capillary blood.


- Adjustable penetration depth.
- Trigger position indicates the before and after use.
- Pre-loaded contact activated safety lancet.
- Quality tri-bevel needle with high-speed penetration.
- Sterile, single use and disposable.
- Ergonomic case and user-friendly design.
- Safety, after use the lancet retracts inside the case, avoiding accidents.
- Easier twist-off cap design.

The Vitrex Flex 3 Adjustable Safety Lancet has three penetration settings:

- 1.3 mm
- 1.8 mm
- 2.3 mm

The dimensions:

- Length: 49.0 ± 1.0 mm
- Width: 17.0 ± 0.3 mm






MATERIALS:

The materials used to manufacture our Vitrex Flex 3 safety lancets are Medical grade:

Component	Material
Needle	Stainless Steel (06Cr19Ni10)
Shell	ABS – QS-0613-016
Trigger	ABS – QS-0613-016

TECHNICAL SPECIFICATIONS:

Colour Code		Configuration	Penetration Depth	Article No.	Application
Blue		23G	1.3 mm 1.8 mm 2.3 mm	395813	Versatile adjustable lancet for most common use.
Yellow		26G	1.3 mm 1.8 mm 2.3 mm	396813	Versatile adjustable lancet for most common use.
Orange		28G	1.3 mm 1.8 mm 2.3 mm	397813	Versatile adjustable lancet for most common use.

Shelf life 5 years.

PRINT:

The terminology, graphic symbols and information applied to this medical device comply with EN 1041:2013 and ISO 15223-1:2016.

PACKAGING:

The inner box is made of white cardboard and the outer carton is made of corrugated cardboard.

All products are packaged according to ISO 11607-1:2019 for medical devices.

Each box contains batch number, expiration date and sterilization information.

- An inner box contains 100 individual Vitrex safety lancet units.
- One outer carton contains 10 inner boxes.



Transit packaging

Pallet

Inner Box

Length	Width	Height	Quantity
13.5 cm	9 cm	9 cm	100 pcs.

Packaging:

Primary packaging:	100 pcs.
Transit packaging:	1,000 pcs.
Pallet:	84,000 pcs.



Vitrex Flex 3

Depth Adjustable Safety Lancet

Technical Data Sheet



LANGUAGES:

The labels contain 13 languages and may be sold in countries where these languages are accepted by national law.

- | | | |
|-----------|--------------|-----------|
| • Danish | • Swedish | • English |
| • Spanish | • German | • French |
| • Italian | • Portuguese | • Russian |
| • Arabic | • Finnish | • Dutch |
| • Greek | | |

QUALITY:

- Quality system in compliance with ISO 9001: 2015 and ISO 13485: 2016.
- Complies with Annex II, (excluding section 4) of the Directive 93/42 EEC.

CHARACTERISTICS:

- Shelf life-time: 5 years.
- Disposable.
- Sterilization by gamma radiation.



Vitrex Medical a/s
Vasekaer 6-8 - DK 2730 Herlev
Phone +45 4494 7011
vitrex@vitrex.dk
www.vitrexmedical.com

