



Forhandler:

MEDI PLAST


Medioplast A/S

Marielundvej 46E, 2730 Herlev

T 43 44 40 00

info.dk@mediplast.com

www.mediplast.com



Technical sheet

YANKAUER Rigid Classic

Reference(s)

Y1200, Y1201, Y1500, Y1501, Y1800, Y1801, Y2200, Y2201, Y2800, Y2801

Description

Classic and rigid suction cannula with:

- Lateral eyes for controlled atraumatic suction
- Beveled tip edges to minimize tissue trauma
- With (code: Yxx01) or without (code: Yxx00) vacuum control to manage the suction
- Tapered universal proximal end for tubing connection

Code	Description	Length (±15 mm)	Diameter (CH)	Holes	Vacuum Control	Pack factor
Y1200	Fine tip (pediatric, ENT use)	275mm	12CH	2	NO	50
Y1201		275mm	12CH	2	YES	50
Y1500	Regular tip (pharyngeal use)	195mm	15CH	0	NO	100
Y1501		195mm	15CH	0	YES	100
Y1800	High tip (orthopedic use, significant particles)	250mm	18CH	6	NO	50
Y1801		250mm	18CH	6	YES	50
Y2200	High tip (orthopedic use, significant particles)	275mm	22CH	6	NO	50
Y2201		275mm	22CH	6	YES	50
Y2800	High tip (orthopedic use, significant particles)	280mm	28CH	4	NO	50
Y2801		280mm	28CH	4	YES	50

Main materials

PVC - DEHP free – latex free

Pack factor/Packaging

Sterile double wrapped individual peel pouch

Sterilization

Sterilized by Ethylene Oxide

Manufacturer

GC MEDICA ENTREPRISE LTD - ISO 13485

Distributed by Medline International France SAS.

Regulation

In compliance with Medical Devices Directive 93/42/EEC

MDD Class: IIa sterile



0197



Latex free



© 2016 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-YANKCL-
EN01

01.01.INT.12