



Laerdal
helping save lives



The BAG™

Single-patient use BVM resuscitator

The BAG™ is a single-patient use, manual resuscitator made with high-quality materials. It is equipped with a range of features that make it userfriendly, such as an integrated hand strap, excellent recoil and a rotatable bag. It is designed to meet the latest safety and performance requirements.

The BAG is available in 3 sizes - Adult, Pediatric and Newborn.

The BAG offers additional features compared to the previous disposable BVM resuscitator, The BAG II.

The BAG™ Product Features

- ✓ Self-inflating bag is made with materials that are both PVC free and Phthalate free*
- ✓ It is designed to meet the latest safety and performance requirements.
- ✓ It has an integrated hand strap that helps to form a secure grip.
- ✓ Highly visible Pressure Relief Valve for ease of use. Adult and Pediatric versions have a simple and intuitive locking system.
- ✓ Made with high quality materials, it offers excellent recoil for a smooth bagging experience.
- ✓ It comes with Laerdal's well trusted face mask, known for its good seal that reduces chance of mask leakage.
- ✓ The BAG can rotate, allowing for uninterrupted ventilation in case you need to change positions.
- ✓ High-quality large valves enable delivery of oxygen concentration of up to 98%.



Integrated PEEP-port and Manometer-port



Clearly visible Pressure Relief Valve on all sizes



Integrated hand strap that helps to form a secure grip.



Accessories: Breathing Filter and Bag Refill Valve Adapter with 25 mm connector.



Accessories: PEEP Valves and Manometer

Ordering Information

The BAG resuscitator (packs of 12)

84711040	The BAG resuscitator Newborn w/mask #1
84722040	The BAG resuscitator Paediatric w/mask #2
84723040	The BAG resuscitator Paediatric w/mask #3
84734040	The BAG resuscitator Adult w/mask #4
84735040	The BAG resuscitator Adult w/mask #5

The BAG Masks. Remove Disposable (packs of 12)

847011	The BAG Mask #1. Remove Disposable
847012	The BAG Mask #2. Remove Disposable
847013	The BAG Mask #3. Remove Disposable
847014	The BAG Mask #4. Remove Disposable
847015	The BAG Mask #5. Remove Disposable

*The oxygen tube and the masks are not PVC free.

Accessories (packs of 12)

847021	The BAG PEEP Valve (2-10 cmH ₂ O). Remove Disposable
847022	The BAG PEEP Valve (5-20 cmH ₂ O). Remove disposable.
847030	The BAG Breathing Filter: Remove disposable.
847040	The BAG Manometer: Remove disposable
847050	The BAG Bag Refill Valve Adapter: Remove disposable.

Regulatory standards and approvals

Europe: CE marked as class IIa device in accordance with Regulation (EU) 2017/745 for medical devices.

Laerdal Medical is certified by DNV Product Assurance AS to ISO 13485:2016

CE2460 The product is in compliance with the essential requirements of Regulation (EU) 2017/745 for medical devices.



For more information, scan the QR code or visit [Laerdal.com](https://www.laerdal.com)

©2025 Laerdal Medical. All rights reserved. 2503-1610204763



Laerdal
helping save lives