




REDIGAUZ[®]

YANGA[®]



TUBULAR
BANDAGE





20 m 100% Cotton



SIZE	∅	CODE		APPLICATION
------	---	------	--	-------------

1,5	1,5 cm	GAUZ-20-1,5-B	GAUZ-20-1,5-N	fingers or toes
2,5	2,5 cm	GAUZ-20-2,5-B	GAUZ-20-2,5-N	bulky fingers or toes dressings
3,5	3,5 cm	GAUZ-20-3,5-B	GAUZ-20-3,5-N	infant arms or fingers bandages with splints
6	6 cm	GAUZ-20-6-B	GAUZ-20-6-N	infant legs, arms or adult legs
8	8 cm	GAUZ-20-8-B	GAUZ-20-8-N	adult legs or infant head
10	10 cm	GAUZ-20-10-B	GAUZ-20-10-N	adult legs or head
12	12 cm	GAUZ-20-12-B	GAUZ-20-12-N	large thighs or adult head
16	16 cm	GAUZ-20-16-B	GAUZ-20-16-N	infant trunk
22	22 cm	GAUZ-20-22-B	GAUZ-20-22-N	adult trunk

⊙ REFERENCE COLOUR	B	N



20 m 100% Syntetic (S)

20 m 100% Polypropylene (Po)



SIZE	∅	CODE		APPLICATION
------	---	------	--	-------------

1,5	1,5 cm	GAUZ-20-1,5-S	GAUZ-20-1,5-Po	fingers or toes
2,5	2,5 cm	GAUZ-20-2,5-S	GAUZ-20-2,5-Po	bulky fingers or toes dressings
3,5	3,5 cm	GAUZ-20-3,5-S	GAUZ-20-3,5-Po	infant arms or fingers bandages with splints
6	6 cm	GAUZ-20-6-S	GAUZ-20-6-Po	infant legs, arms or adult legs
8	8 cm	GAUZ-20-8-S	GAUZ-20-8-Po	adult legs or infant head
10	10 cm	GAUZ-20-10-S	GAUZ-20-10-Po	adult legs or head
12	12 cm	GAUZ-20-12-S	GAUZ-20-12-Po	large thighs or adult head
16	16 cm	GAUZ-20-16-S	GAUZ-20-16-Po	infant trunk
22	22 cm	GAUZ-20-22-S	GAUZ-20-22-Po	adult trunk

S	Po
SYNTHETIC	POLYPROPYLENE

Redigauz tubular bandage represents an extremely practical and reliable method of fixing wound dressing.

Tubular bandages replace the need for traditionally used flat gauze bandages which must be wound around the part to be covered and permit full freedom of movement.


Redigauz bandages have an unbroken surface, and therefore do not come loose; wounds or sensitive areas of skin remain protected.

Redigauz is highly versatile, being suitable also for traction bandages, as a skin shield under plaster of Paris bandages, or as a covering material for splints and cotton wool padding.

Redigauz is manufactured from bleached cotton, natural cotton, synthetic and polypropylene.

It's water repellent.

Redigauz can be autoclaved to provide sterile dressings.



The knitted fabric.



Redigauz is water repellent.



REDIGAUZ[®]

Ø 8



Ø 10



SANITARY DIVISION OF
YANGA[®]
ITALIAN TEXTILE GROUP



MANUFACTURED BY YANGA
ITALY



COMPANY
WITH QUALITY SYSTEM
CERTIFIED BY DNV
= ISO 9001/2000 =

YANGA SRL		13864 CREVACUORE [BIELLA] ITALY		TEL +39 015 768256		FAX +39 015 768910
		E-MAIL yanga@yanga.it		WEBSITE www.yanga.it		