


PRODUCT DATA SHEET

Tungespartel russel-davis 85x38mm (OM173R)

RUSSEL-DAVIS Tunge Depressor, klinge, m/ rille til endotrach. tubus, 38 x 85 mm, ikke-steril, genanvendelig



Attribute Name	Attribute Value
Product Search Type	Mouth Gag
Author	RUSSEL-DAVIS
Materiale	General Surgical Instruments, ENT
Type	Tongue Depressor
Type	blade only
Længde	135.00 mm
Length (Inch)	5 1/4"
Blades / Jaw	w/groove for endotrach.tubes
Jaw/Blade width	38.00 mm
To be used with (code)	OM104R, OM119R, OM105R, OM116R, OM118R, OM169R
Facon / form	Ear-Nose-Throat Surgery, Maxillo-Facial Surgery
Facon / form	Tonsillectomy, Adenoidectomy
UDI direkte mærkning	Yes

Attribute Name	Attribute Value
Materiale	reusable
Type	non-sterile
Type	No
Jaw/Blade depth	85.00 mm