

Careline

Leg bags, Sterile



Intended use:

Urine drainage leg bags with measuring scale for collection and measuring of urine output, fitted with bottom outlet.

List of product range:

Volume	Tube Length	KombiKon	Sample Port	Cross Valve	Lever Tap	Hospital Pack	Community Pack	REF
500ml	10cm	Yes	-	Yes	-	Yes	-	12 618 1S
500ml	10cm	Yes	-	-	Yes	Yes	-	12 629 1S
500ml	27cm	Yes	-	Yes	-	Yes	-	12 623 1S
500ml	27cm	Yes	-	-	Yes	Yes	-	12 632 1S
350ml	10cm	-	Yes	Yes	-	-	Yes	08 613 1S
350ml	27cm	-	Yes	Yes	-	-	Yes	08 620 1S
500ml	10cm	-	Yes	Yes	-	-	Yes	12 613 1S
500ml	10cm	Yes	-	-	Yes	-	Yes	12 612 1S
500ml	27cm	Yes	-	Yes	-	-	Yes	12 620 1S
500ml	27cm	Yes	-	-	Yes	-	Yes	12 621 1S
750ml	10cm	Yes	-	Yes	-	-	Yes	24 610 1S
750ml	10cm	Yes	-	-	Yes	-	Yes	24 612 1S
750ml	27cm	Yes	-	Yes	-	-	Yes	24 620 1S
750ml	27cm	Yes	-	-	Yes	-	Yes	24 621 1S

Environmental issues:

Contains Latex?.....No
 Contains Paper?.....Yes
 Contains Cotton?.....No
 Tube inlet/outlet, Foil.....Polyvinyl Chloride
 KombiKon/Sample Port...Polyvinyl Chloride, Polypropylene, Styrene Butadiene Styrene
 Hanger, Tube Clamp.....Polypropylene
 Cross Valve.....Polyethylene, Polypropylene
 Cap.....Polyethylene
 Silicone Tube.....Silicone
 Lever Tap.....Polypropylene

Careline Leg Bags, Sterile

Standardisation:

CE-marking. **CE 0459**
 Quality system certified to **ISO 13485**
 GMDN. **31085**

Sterilisation:

Sterile. **Yes**
 Sterilised by **Ethylene Oxide**
 Shelflife **5 years from the date of sterilisation**

Instruction for use:

Yes ☒ ■
 No ☐ ■

Storage and handling:



Packaging:

REF	Shipping carton				MOQ, pcs
	Box Qty, pcs	Dimensions, mm	Gross Wt, kg	Volume, m ³	
08 613 1S	100	792x218x392	6,6	0,067	2800
08 620 1S	100	792x218x392	6,6	0,067	2800
12 612 1S	100	792x218x392	6,0	0,067	2000
12 613 1S	100	792x218x392	6,3	0,067	2000
12 618 1S	100	792x218x392	5,6	0,067	2000
12 620 1S	100	792x218x392	7,7	0,067	2000
12 621 1S	100	792x218x392	7,9	0,067	2000
12 623 1S	100	792x218x392	6,4	0,067	2000
12 629 1S	100	792x218x392	6,0	0,067	2000
24 610 1S	100	792x218x392	7,1	0,067	2000
24 612 1S	100	792x218x392	7,1	0,067	2000
24 620 1S	100	792x218x392	6,9	0,067	2000
24 621 1S	100	792x218x392	7,1	0,067	2000