


DANSU A/S

PRODUCT SPECIFICATION SHEET FOR PRODUCT NO. 50-200

Description : Wood hospital applicator, 150 x 2.2mm
Cotton head, Bag 100 pcs
10 bags/box – 20 boxes/carton – 24 cartons/pallet



Applicator:

measurements: A: 150.0 ± 2.0 mm

Stick:

measurements: B: 146.0 ± 1.0 mm

C: 2.2 ± 0.1 mm

material: White birchwood, completely round, straight and polished

Head:

measurements: D: 5.5 ± 0.5 mm

E: 14.0 ± 1.0 mm

material: Cotton – pharmaceutical grade

weight: 0.04 gr ± 10%

Adhesive:

glue: Gludan non-toxic

Packaging:

bag material: Polyethylene PE-LD 04
size: 110x210x0.04 mm, with holes

carton: material: Corrugated cardboard, no. 29
size: 365x350x460 mm
gross weight: 10.7 kg