


71-504

TECHNICAL DATA FOR : PRODUCT NUMBER 143-2504

Description : PP applicator 400 x 4 mm - non sterile
Mega cotton head



- Indication** : Specimen collection
: Single use
- Applicator** : measurement A : 405 ± 1 mm
- Stick** : measurement B : 400 ± 1 mm
: measurement C : 4 ± 0.1 mm
: material : White Polypropylene
- Head** : measurement D : 14 ± 1 mm
: measurement E : $35 \text{ mm} \pm 2$ mm
: material : 100% Cotton (European Pharmacopea)
- Fixing method** : Heat melting system + non toxic glue
- Packaging** : 50 applicators/bag
: 10 bags/carton
: 120 cartons/pallet
- Carton n° 7** : size : 500 x 210 x 135 mm
: gross weight : 2.3 kg
- Handling and storage** : normal conditions
- Environment** : harmless
- Labelling** : **CE**

Manufactured in Switzerland by Applimed SA. Iso 9002/EN46002	Edition doc.: 02.98 File : c\wp51\143-2504 Révision : 1	Visa > Date >	Issued by Vérfifié le 25 MAI 1998	Controlled by
--	---	------------------	---	---------------

This drawing contains confidential and proprietary information. No use or disclosure, either in whole or in part, may be made without authorization by APPLIMED SA.