


Product Data Sheet:

Spikes

CODAN Spike/SWAN



1. A compact and user friendly system for preparing and withdrawing medicines. Specially shaped piercing device. The pressure-resistant hydrophobic 0.2µm filter guarantees sterile ventilation and reliably excludes toxic aerosols during the preparation of medicines. The high-quality housing materials ensure that there is no loss of active substances during the preparation and withdrawal of medicines. The large surface area of the filter membrane ensures easy fluid transfer without high working pressure. Furthermore, fluids can be withdrawn without leaving significant residues.
2. SWAN = Swabable adapter Needle free adapter with Luer activated open/close mechanism (Valve), with female Luer-Lock connectors. For giving or taking of solutions needle free.

Ref.	Product Description	Packaging Units
16.5290	CODAN Spike + needlefree adapter	50/600/14400
16.5287	25x CODAN Spike + needlefree adapter	25/500/31500