

care*fix*[®]
Creating functional comfort



HEAD fixation

Precut elastic dressing fixation



TYTEX[®]

The CareFix head fixations have been designed to hold bandages or surgical dressings in place without tape. They are precut, easy to use and time-saving for caregivers.

Features:

- Effective retention of all types of dressings
- Precut and ready to use
- Simple and fast application
- Elasticity and support
- Comfortable fit
- Fit children and adults
- Available in 4 flexible sizes
- Practical colour-coding indicates sizes
- Washable (60°C/140°F)
- Latex-free material for the sake of the patient and caregiver
- Tested for harmful substances according to Oeko-Tex standard 100

Indications:

For usage after treatment of ear, face, head, and neck disorders

Composition:

Head: 60% Polyamide, 18% Elastane, 16% Cotton, 6% Polyester.

Headband: 50% Cotton, 22% Elastane, 18% Polyamide, 10% Polyester.

Washing instructions:



Recommended sizes:

Head:

S = Child, M = Child/adult, L = Adult, XL = Adult

Headband:

One Size

Packaging:



Dispenser boxes w/10 pcs.



Bag w/10 pcs.



For further information please contact:

Dealer: