


carefix[®]

Creating functional comfort



TUBE fixation

Precut IV and dressing fixation



The CareFix tubes are used for application to the arm and leg, for IV fixation and for retaining surgical dressings and bandages.

The tubes are precut, which makes them ideal for emergency rooms and hospital units where time is essential.

Features:

- Hold bandages or surgical dressings in place without tape
- Quick and easy to apply
- Flexible, breathable, comfortable and soft to the skin
- Timesaving, ready-to-use precut lengths
- Designed to provide easy access to the injection valve on catheters
- With precut opening for the thumb in sizes S and M
- Can be cut according to need without the material running
- Practical colour-coding indicates sizes
- Washable
- Latex-free material for the sake of the patient and caregiver
- Tested for harmful substances according to Oeko-Tex standard 100

Indications:

- IV fixation
- Dressing fixation
- Bandage fixation

Composition:

70% Polyamide, 21% Elastane, 9% Polyester.

Washing instructions:



60°C/140°F

Recommended sizes:

S = Arm (child)

M = Arm (child/adult)

L = Arm/leg (adult)

XL = Leg (adult)

Packaging:



Dispenser boxes w/15 pcs.



For further information please contact:

Dealer:



TYTEX A/S has been certified according to ISO 9001, ISO 14001 and ISO 13485

