

Technical Data Sheet

Sorbact® Ribbon Gauze



General	
Product description	A green cotton fabric with permanently bound dialkyl carbamoyl chloride (DACC)
Ref no	Size
98118	1x50 cm
98119	10x200 cm
98120	5x200 cm
98121	2x50 cm

Material		
Layer 1	Bacteria binding layer	A green cotton fabric with permanently bound dialkyl carbamoyl chloride (DACC)

Packaging	
Primary packaging	Individual hot sealed paper/paper pouch
Secondary packaging	Cardboard box
Tertiary packaging	Well transport box

Manufacturing and regulatory	
Sterilization method	Gamma irradiation
Shelf life	5 years
Storage conditions	Keep dry and away from direct sunlight.
Medical Device Class acc to Annex IX of 93/42/EEC	IIb
Country of manufacture	Sweden
Certification	ISO 13485:2012 ISO 14001:2004 ISO 9001:2008

Logistics				
Ref no	98118	98119	98120	98121
Weight, unit (kg) (device in primary packaging)	0.002	0.011	0.006	0.002
Weight, secondary packaging (kg) (with device)	0.086	0.163	0.110	0.089
Weight, tertiary packaging (kg) (with device)	0.8	1.2	0.9	0.8
Units in secondary packaging (pcs)	20	10	10	20
Units in tertiary packaging (pcs)	120	60	60	120
EAN code primary packaging	7392130181180	7392130181197	7392130181203	7392130181210
EAN code secondary packaging	7392130981186	7392130981193	7392130981209	7392130981216
EAN code tertiary packaging	7392130800333	7392130800326	7392130800074	7392130800050

Product specifics	
Package material: bleach method (chlorine, ECF, TCF)	Non chlorine
Package material: PVC or other halogen based plastic material	No
Addition of antimicrobial substance	No
Addition of antibiotics	No
Addition of allergens:	
-Latex	No
-Lanolin	No
-Colophony and derivatives	No
-Phthalates (DEHP, DINP, DIDP, BBP, DBP)	No
-Triclosan	No