

Technical Data Sheet

Sorbact[®] Gel Dressing



General	
Product description	Bacteria binding gel dressing: acetate fabric with permanently bound dialkyl carbamoyl chloride (DACC) and gel containing carbomer, sodium hydroxide, propylene glycol and demineralized water
Ref no	Size
98136	7.5x7.5 cm
98137	7.5x15 cm

Material		
Layer 1	Bacteria binding layer	A green acetate fabric with permanently bound dialkyl carbamoyl chloride (DACC)
Layer 2	Gel	Carbomer Propylene glycol Sodium hydroxide Demineralized water

Packaging	
Primary packaging	Individual hot sealed alu/alu foil
Secondary packaging	Cardboard box
Tertiary packaging	Well transport box

Manufacturing and regulatory	
Sterilization method	Beta irradiation
Shelf life	3 years
Storage conditions	Keep dry and away from direct sunlight.
Medical Device Class acc to Annex IX of 93/42/EEC	IIb
Country of manufacture	Sweden
Certification	ISO 13485:2012 ISO 14001:2004 ISO 9001:2008

Logistics		
Ref no	98136	98137
Weight, unit (kg, device in primary pack)	0.009	0.017
Weight, secondary pack (kg, with device)	0.101	0.185
Weight, tertiary pack (kg, with device)	4.4	4.5
Units in secondary pack (pcs)	10	10
Units in tertiary pack (pcs)	360	210
EAN code primary pack	7392130181364	7392130181371
EAN code secondary pack	7392130981360	7392130981377
EAN code tertiary pack	7392130800234	7392130800241

Product specifics	
Package material: bleach method (chlorine, ECF, TCF)	Non chlorine
Package material: PVC or other halogen based plastic material	No
Addition of antimicrobial substance	No
Addition of antibiotics	No
Addition of allergens:	
-Latex	No
-Lanolin	No
-Colophony and derivatives	No
-Phthalates (DEHP, DINP, DIDP, BBP, DBP)	No
-Triclosan	No