

Technical Data Sheet

Sorbact[®] Compress



General	
Product description	Bacteria binding compress (8 layers): acetate fabric with permanently bound dialkyl carbamoyl chloride (DACC)
Ref no	Size
98128	4x6 cm
98125	7x9 cm

Material		
Layer 1	Bacteria binding layer	A green acetate fabric with permanently bound dialkyl carbamoyl chloride (DACC)

Packaging	
Primary packaging	Individual hot sealed paper/plastic pouch
Secondary packaging	Cardboard box
Tertiary packaging	Well transport box

Manufacturing and regulatory	
Sterilization method	Gamma irradiation
Shelf life	5 years
Storage conditions	Keep dry and away from direct sunlight
Medical Device Class acc to Annex IX of 93/42/EEC	IIb
Country of manufacture	Sweden
Certification	ISO 13485:2012 ISO 14001:2004 ISO 9001:2008

Logistics		
Ref no	98128	98125
Weight, unit (kg) (device in primary packaging)	0.003	0.005
Weight, secondary packaging (kg) (with device)	0.224	0.224
Weight, tertiary packaging (kg) (with device)	1.2	1.6
Units in secondary packaging (pcs)	40	40
Units in tertiary packaging (pcs)	240	240
EAN code primary packaging	7392130181289	7392130181258
EAN code secondary packaging	7392130981285	7392130981254
EAN code tertiary packaging	7392130800029	7392130800043

Product specifics	
Package material: bleach method (chlorine, ECF, TCF)	Unbleached, ECF or TCF
Package material: PVC or other halogen based plastic material	Absent
Addition of antimicrobial substance	No
Addition of antibiotics	No
Addition of allergens:	
-Latex	No
-Lanolin	No
-Colophony and derivatives	No
-Phthalates (DEHP, DINP, DIDP, BBP, DBP)	No
-Triclosan	No