

Technical Data Sheet

Sorbact[®] Foam Dressing



General	
Product description	Acetate fabric with permanently bound dialkyl carbamoyl chloride (DACC) laminated on polyurethane foam
Ref no	Size
98310	10x10 cm
98315	15x15 cm
98320	10x20 cm


Material		
Layer 1	Foam	Absorbent hydrophilic polyurethane foam 5 mm with polyurethane backing film
Layer 2	Adhesive	Copolyester
Layer 3	Bacteria binding layer	A green acetate fabric with permanently bound dialkyl carbamoyl chloride (DACC)

Packaging	
Primary packaging	Individual hot sealed paper/paper pouch
Secondary packaging	Cardboard box
Tertiary packaging	Well transport box

Manufacturing and regulatory	
Sterilization method	Gamma irradiation
Shelf life	3 years
Storage conditions	Keep dry and away from direct sunlight
Medical Device Class acc to Annex IX of 93/42/EEC	IIb
Country of manufacture	Sweden
Certification	ISO 13485:2012 ISO 14001:2004 ISO 9001:2008

Logistics			
Ref no	98310	98315	98320
Weight, unit (kg, device in primary pack)	0.019	0.035	0.032
Weight, secondary pack (kg, with device)	0.16	0.30	0.27
Weight, tertiary pack (kg, with device)	1.5	2.8	2.6
Units in secondary pack (pcs)	10	10	10
Units in tertiary pack (pcs)	80	80	80
EAN code primary pack	7392130183108	7392130183153	7392130183207
EAN code secondary pack	7392130983104	7392130983159	7392130983203
EAN code tertiary pack	7392130800364	7392130800395	7392130800388

Product specifics	
Package material: bleach method (chlorine, ECF, TCF)	Unbleached, ECF or TCF
Package material: PVC or other halogen based plastic material	No
Addition of antimicrobial substance	No
Addition of antibiotics	No
Addition of allergens:	
-Latex	No
-Lanolin	No
-Colophony and derivatives	No
-Phthalates (DEHP, DINP, DIDP, BBP, DBP)	No
-Triclosan	No



Layer 1: Foam with backing film

Layer 2: Adhesive

Layer 3: Bacteria binding layer