



A multi-purpose couch providing electric positioning for gynaecology procedures, as well as a full length couch surface for most general prone and supine examinations and minor treatment procedures

## MEDICARE GP GYNAECOLOGY COUCH

### MEDICARE GP GYNAECOLOGY COUCH

**SM8543** - Electric height, electric tilt, electric backrest and manual head section



Model **SM8543** with Lithotomy Stirrups. Colour Black

### STANDARD MODEL SPECIFICATION

- Safe working load and lifting capacity (from minimum height) of 240kgs (528lbs=37.7 stone)
- Electric adjustable angle backrest from horizontal to +80°
- Electrically adjustable seat angle with trendelenburg facility to -10°
- Power assisted adjustable head section from -25° to +50°
- Breathing hole and plug fitted as standard
- Excellent height range from 50-103cm (20"-41")
- Retractable wheel design base frame with adjustable foot for uneven floor
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard
- Supplied complete with Lithotomy stirrups and debris tray
- Standard width 70cm (28")
- Paper towel holder fitted as standard

### OPTIONAL COUCH ACCESSORIES

- Battery backup
- Gynae Leg Supports
- Head Support Cushion
- Triple Foot Switch
- Matching operators chair or stool



Model **SM8543** with raised head section



Model **SM8543** with head section negative



Model **SM8543** with electric tilted seat angle



Model **SM8543** at Low transfer height