

A very strong and rigid couch available with a wide range of optional accessories. Suitable for supporting very heavy patients for examination and minor treatment procedures.



2 SECTION MEDICARE COUCH

2 SECTION MEDICARE COUCH CHOICE OF MODELS AVAILABLE


SM2550 Hydraulic

SM2560 Electric

SM2570 Electric with electric backrest



Width 65cm



Model **SM2560**. Colour Cream

STANDARD MODEL SPECIFICATION

- Safe working load and lifting capacity (from minimum height) of 240kgs (528lbs=37.7 stone)
- Adjustable backrest from -25 to +75 degrees positive (0 to +85 degrees Model SM2570)
- The power assisted backrest section incorporates anti-trap design (model SM2550/2560)
- Choice of manual or electric height elevation
- Excellent height range from 43-97cm (17"-37")
- Retractable wheel design and adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Standard width 65cm (25")
- Choice of upholstery colours - see back cover
For alternative widths, please refer to price list.

OPTIONAL COUCH ACCESSORIES

- Model 6053 Protective base frame shrouding
- Model 6004 Breathing hole and plug for prone treatments
- Model 6200 Single foot switch operation for model SM2560
- Model 6207 Double foot switch operation for model SM2570
- Model 6005 Paper towel holder
- Model 6002 Side support rails (each)
- Model 6002P Side support rail cushions (each)
- Model 6006 Upper limb support
- Model 6007 Removable head/knee support
- Model 6075 Optional 75cm width
- Model 6085 Optional 85cm width

See pages 34 and 35 for additional information



Model **SM2560**
Showing negative backrest angle.
Colour Cream



Model **SM2560**
Minimum height. Colour Cream