

A simple 2 section couch, with an adjustable backrest for supporting the patient in a sitting position whilst being treated or examined. An optional breathing hole and plug can also be added for prone treatments


MODEL SPECIFICATION

- Safe working load and lifting capacity (from minimum height) of 200kgs (440lbs=31.5 stone)
- Power assisted backrest angle adjustment from -25 degrees to +75 degrees
- Excellent height range from 46-99cm (18"-39")
- 100mm diameter individual locking castors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard
- Standard width 64cm (25")
- Choice of upholstery colours - Dark Grey, Dark Blue or Black

OPTIONAL COUCH ACCESSORIES

- Breathing hole and plug
- Single foot switch operation
- Paper towel holder
- Matching therapy chair or stool
- Removable head/knee support
- Stretch towelling covers

STERLING 2 SECTION COUCH - Model No. SX2007



Model **SX2007** - with optional breathing hole. Head section in negative position. Colour Dark Blue



Model **SX2007** - with optional breathing hole. Low transfer height. Colour Dark Blue



Model **SX2007** - with optional breathing hole. Colour Dark Blue.