

CODAN Adapter

SWAN-LOCK®



1. **SWAN = SWabater Needle-free Adapter** with Luer activated open/close mechanism (Valve), with Male and Female Luer-Lock connectors. For needle-free administration or taking of solutions. Variant - **SWAN-LOCK® RED** with red transparent housing, colour coded for identification of e.g. arterial access.
Technical Data:
 - Priming volume: 0,09 ml
 - Negative displacement: 0,02 ml
 - Flow reduction: 16%
 - Bi-directional flow, for addition/withdrawal
 - Pressure resistance: 200 kPa according to EN ISO 8536-8
 - Negative pressure resistance: -20 kPa according to EN ISO 8536-8
 - No protection cap is required for the Female Luer-Lock connector, capless version
 - Transparent housing
 - Free of Latex and PVC
 - Resistant to lipids/alcohol
 - Resistant to cytotoxic agents tested with: Vespesid/Bristol Arzneimittel GmbH and Endoxan/Asta Medica AG
 - MRI compatible
 - Microbial barrier tested, 140 activations over a 7-day period



Disclaimer:

The above is indicated only for information. This document is not a specification. The values given are subject to the usual fluctuations. Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times.

Subject to alterations without prior notice!

REF	Product description	Packaging	GTIN	Filling/Storage volume
16.5267	Swabable needle-free adapter SWAN-LOCK®	1 100 3200 76800	04055698002511 04055698002528 04055698002535 -	



CODAN Adapter

SWAN-LOCK®

REF	Product description	Packaging	GTIN	Filling/Storage volume
16.5273	Swabable needle-free adapter SWAN-LOCK® RED	1	04055698002542	
		100	04055698002559	
		3200	04055698002566	
		76800	-	