

Technical Description

Connecting Tube for Drainage Bag

| | |
|-------------------------------------|---|
| Trade Name | Connecting Tube for Drainage Bag |
| Manufacturer | PBN MEDICALS Denmark A/S Knud Bro Alle 3 DK-3660 Stenlose, DENMARK |
| Product Description | Connecting tubes for all types of drainage bags. Available with or without stopcocks. |
| Classification | Is |
| GMDN Code | 45389 |
| Patent Number (if any) | None |
| EC Certificate | DGM-375 |
| ISO Certificate | DGM-377 |
| Notified Body | DGM (0543) |
| Contains | Connecting tube |
| Material | Hostaform® C 9021 nature and PVC |
| Facility Certification | PBN MEDICALS is FDA registred and certified in accordance with ISO 13485, the Medical Device Directive 93/42/EEC and the canadian ISO 13485. |
| Quality Control | Every lot of components is inspected before entering production. The finish products are subject to 100% visual control. |
| Method of Sterilization | Ethylene Oxide |
| Sterilization Facility | Sterigenics Germany GmbH, Rheingaustraße 190-196, 65203 Wiesbaden, Germany |
| Packaging | <u>Protected Unit:</u> Individually sealed in uncoated Perfecseal® pouch and packed in a cardboard box <u>Quantity per box:</u> Box of 25 <u>Label:</u> Conform to EN 980 and EN 1041 standards |
| Shelf life | 5 years |
| Recommendations | Sterility guaranteed for undamaged and unopened packaging prior to use Use before the expiration date on labels Store in a dark, dry, cool place; Avoid extended exposure to light |
| Information related to Latex | Free of Latex |