

Oxygen Masks

High Concentration

Intended use:

High concentration masks intended for administration of oxygen. They are equipped with a reservoir bag allowing an adequate amount of oxygen to be available to meet unpredictable breathing patterns and tidal volumes. Oxygen concentration delivered is 80%-100% depending on the patient's breathing rate and tidal volume. The product is intended for single use only.



List of product range:

Oxygen mask, High concentration, Elongated (under the chin) with adjustable nose-clip and elastic strap				Oxygen tubing		
	REF*	Size	Description	m	Type	Connector
	104MM	Adult	Rebreathing mask, Reservoir bag, with tubing	2.1	Standard lumen	Standard
	9104MM	Adult	Rebreathing mask, Reservoir bag, with tubing	2.1	No-crush lumen	Standard
	105MM	Adult	Rebreathing mask, Reservoir bag, no tubing	-	-	-
	106MM	Adult	Non-rebreathing mask, Reservoir bag, Check valve, One side valve, with tubing	2.1	Standard lumen	Standard
	107MM	Adult	Non-rebreathing mask, Reservoir bag, Check valve, One side valve, no tubing	-	-	-
	108MM	Adult	Non-rebreathing mask, Reservoir bag, Check valve, Two side valves, with tubing	2.1	Standard lumen	Standard
	9108MM	Adult	Non-rebreathing mask, Reservoir bag, Check valve, Two side valves, with tubing	2.1	No-crush lumen	Standard
	109MM	Adult	Non-rebreathing mask, Reservoir bag, Check valve, Two side valves, no tubing	-	-	-
	3224MM	Pediatric	Rebreathing mask, Reservoir bag, with tubing	2.1	Standard lumen	Fits-all
	225MM	Pediatric	Rebreathing mask, Reservoir bag, no tubing	-	-	-
	226MM	Pediatric	Non-rebreathing mask, Reservoir bag, Check valve, One side valve, with tubing	2.1	Standard lumen	Standard
	3226MM	Pediatric	Non-rebreathing mask, Reservoir bag, Check valve, One side valve, with tubing	2.1	Standard lumen	Fits-all
	3230MM	Pediatric	Non-rebreathing mask, Reservoir bag, Check valve, Two side valves, with tubing	2.1	Standard lumen	Fits-all
	93230MM	Pediatric	Non-rebreathing mask, Reservoir bag, Check valve, Two side valves, with tubing	2.1	No-crush lumen	Fits-all



Oxygen mask, High concentration, Standard (under the lip) with adjustable nose-clip and elastic strap				Oxygen tubing		
	REF*	Size	Description	m	Type	Connector
	204MM	Adult	Rebreathing mask, Reservoir bag, with tubing	2.1	Standard lumen	Standard
	205MM	Adult	Rebreathing mask, Reservoir bag, no tubing	-	-	-
	3264MM	Infant	Rebreathing mask, Reservoir bag, with tubing	2.1	Standard lumen	Fits-all
	265MM	Infant	Rebreathing mask, Reservoir bag, no tubing	-	-	-
	3266MM	Infant	Non-rebreathing mask, Reservoir bag, Check valve, Two side valves, with tubing	2.1	Standard lumen	Fits-all
	93266MM	Infant	Non-rebreathing mask, Reservoir bag, Check valve, Two side valves, with tubing	2.1	No-crush lumen	Fits-all
	267MM	Infant	Non-rebreathing mask, Reservoir bag, Check valve, Two side valves, no tubing	-	-	-



* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex?	No
Contains DEHP?	Yes
Mask	Polyvinyl chloride
Elastic strap	(Latex-Free) Polyester
Nose clip	Aluminum
Mask connector/Nipple	High density polyethylene
Reservoir bag	Polyvinyl chloride
Check valve	Oxygen "T" Connector: Polypropylene Flapper: Silicone (non-latex)
One way valve	Flapper: Silicone (non-latex)
Oxygen tubing	Polyvinyl chloride
Tube connector	Polyvinyl chloride

Certification:

CE-marking?	Yes (CE 0120)
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485
Global Medical Device	
Nomenclature (GMDN)	35171/35173

Sterilisation:

Sterile? No
Shelf life 5 years

Oxygen Masks

High Concentration

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
104MM	*	*	50	355 x 406 x 305	4.54	0.044	1,250
105MM	*	*	50	406 x 355 x 254	3.63	0.037	1,250
106MM	*	*	50	355 x 406 x 305	4.76	0.044	1,250
107MM	*	*	50	381 x 254 x 203	3.63	0.029	1,250
108MM	*	*	50	355 x 406 x 305	4.76	0.044	1,250
109MM	*	*	50	406 x 355 x 254	2.49	0.037	1,250
204MM	*	*	50	355 x 406 x 305	4.54	0.044	1,250
205MM	*	*	50	406 x 355 x 254	4.54	0.037	1,250
225MM	*	*	50	406 x 355 x 254	4.31	0.037	1,250
226MM	*	*	50	355 x 406 x 305	4.31	0.044	1,250
265MM	*	*	50	381 x 254 x 203	1.04	0.029	1,250
267MM	*	*	50	381 x 254 x 203	1.81	0.029	1,250
3224MM	*	*	50	355 x 406 x 305	4.31	0.044	1,250
3226MM	*	*	50	406 x 355 x 254	4.54	0.037	1,250
3230MM	*	*	50	355 x 406 x 305	4.54	0.044	1,250
3264MM	*	*	50	406 x 355 x 254	3.40	0.037	1,250
3266MM	*	*	50	406 x 355 x 254	3.40	0.037	1,250
9104MM	*	*	50	355 x 406 x 305	4.54	0.044	1,250
9108MM	*	*	50	355 x 406 x 305	4.99	0.044	1,250
93230MM	*	*	50	355 x 406 x 305	4.76	0.044	1,250
93266MM	*	*	50	406 x 355 x 254	3.86	0.037	1,250

* = Single packed - no shelf box

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Inc.**
5701-1S, Ware Road
Address: McAllen, Texas 78503
USA

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

