


KleenSpec® 590 Series Vaginal Speculum Specifications



Size	590XS - Extra Small		59000 - Small		59001 - Medium		59004 - Large	
*	INCH	CM	INCH	CM	INCH	CM	INCH	CM
A	6.04	15.34	6.94	17.63	6.94	17.63	7.65	19.43
B	0.65	1.63	0.88	2.24	1.16	2.94	1.4	3.56
C	5.48	13.92	5.43	13.79	5.43	13.79	5.45	13.84
D	0.58	1.47	0.74	1.88	0.86	2.18	0.86	2.18
E1	1.95	4.95	2.18	5.54	2.26	5.74	2.52	6.4
E2	2.35	5.97	2.73	6.95	2.89	7.34	3.15	8
F	2.38	3.05	2.68	6.81	2.68	6.81	3.5	8.89
R	0.36	0.91	0.625	1.59	0.625	1.59	0.84	2.13