

SafeTway[®] First

Reduce risk from cross-infection and protect both patient and equipment from contamination with the use of disposable filters & SafeTway[®] mouthpieces

SafeTway[®]

The low cost, low resistance SafeTway[®] mouthpiece allows entry of exhaled air into the measuring system, but not back to the patient. This simple yet effective patented design prevents inhalation of aerosolised or re-aerosolised droplets from subjects previously tested. Fits pulmonary function test equipment, including Peak Flow Meters.

Where biting a mouthpiece of 30mm diameter is difficult, use the 22mm Mini SafeTway[®].



SafeTway[®] mouthpiece



Mini SafeTway[®] mouthpiece

Single-Use



Disposable

BVF

The Bacterial/Viral Filter (BVF) from Vitalograph is specially designed for use on respiratory function testing equipment. It combines a very efficient (99.99%) filter medium (tested at 1 micron) with a very low resistance to flow to ensure that dynamic back pressure remains within ATS standards. For the very old or very young we offer the Mini BVF.



Bacterial/Viral Filter



Mini Bacterial/Viral Filter

Adaptors available to fit most respiratory test equipment




.... from **Vitalograph[®]**
the spirometry people

Mini Mouthpieces



It is important for the subject to lightly bite the mouthpiece in order to obtain a good lip seal. In cases of the very young and particularly the very old, or in cases of facial trauma, a smaller mouthpiece may be required. The Mini Mouthpiece has a 22mm diameter. The adaptor converts from the standard 30mm diameter.



Standard Cardboard Mouthpieces



The traditional single use disposable cardboard mouthpiece is effective in reducing any risk of cross-contamination between subjects by limiting mucosal contact, particularly saliva from the mouth.

The Vitalograph mouthpieces have special linings, inside to prevent cardboard dust inhalation and outside to prevent lip sticking and consequent bleeding.

For your promotional purposes, Vitalograph offer the addition of your logo and a range of pack sizes (normally 50,000 mouthpieces minimum). This is also available on the SafeTway.



Ordering Information

28200Bacterial / Viral Filter (50)
28202Mini Bacterial Viral Filter (30)
20242SafeTway Mouthpieces (Box of 200)
20980Mini SafeTway Mouthpieces (Box of 50)
20201Disposable Cardboard Mouthpieces (200)
20202Mini Mouthpieces (Box of 400)

20302Mini Mouthpiece Adaptor
20303Noseclips (10)
AdaptorsA range to fit most PFT equipment
Please Call.....Customised Mouthpiece Packs, your logo on mouthpiece and/or carton. Variety of pack sizes from 20 m/ps upwards. (50,000 minimum)

Adaptors for all Spirometers/PFMs



The SafeTway, BVF and high quality cardboard mouthpieces from Vitalograph are widely used on Peak Flow Meters, Spirometers and Lung Function Analysers which may have different sizes of connector. Vitalograph has developed a range of adaptors to fit most brands of such devices. If you have an unusual device, no problem! We will custom manufacture one for you at no extra cost. Please call with exact sizes, manufacturer's brand/model or send a sample.



.... from **Vitalograph®**
the spirometry people