



LABORATORIES, INC.

286 ELDRIDGE ROAD, FAIRFIELD, NEW JERSEY 07004 U.S.A.

TECHNICAL DATA SHEET

SPECTRA® 360 ELECTRODE GEL

Description:	Viscous, clear, blue-green, aqueous gel
Composition:	Reverse osmosis (RO) water, humectant, polymers, preservatives and FD&C color
Preservatives:	Propyl paraben and methyl paraben in bacteriostatic concentration
pH Range:	6.3 – 7.0
Viscosity:	175,000 – 260,000 CPS (Brookfield Model RVT, T-C Spindle, 2.5 RPM)
Biocompatibility:	Meets requirements for ISO 10993-1
Precautions:	No special handling required. It is non-flammable, non-corrosive, and is not affected by variations in temperature. Does not contain alcohol
Shelf Life:	Four years from month of manufacture

10-27-10