

# NEX IODIO P2<sup>®</sup> SURGICAL SCRUB

Hand Hygiene is the basis of a correct and effective infection prevention in the hospital environment.

Surgical hand antiseptics, to destroy transient micro-organisms and inhibit the growth of resident micro-organisms, is routinely carried out before undertaking invasive procedures.

Antisepsis may reduce the risk of surgical site infections in patients.



## NEX IODIO P2 7,5% (0,75% free iodine)

### Description

NEX IODIO P2 is a single use brush/sponge with nailcleaner impregnated with Povidone Iodine antiseptic solution.

Expiry: 3 years

### Formulation Composition

Ingredients for 100g.:

Povidone Iodine g. 7,5

Coformulants and water q.s. to g. 100

### Indications

Surgical scrubbing and washing of hands for antiseptic pre-op procedures.

### Microbiological properties

Bactericidal: EN 1040, EN 12054 (Surgical handrub and handwash) and EN 12791 (Surgical hand disinfection)

### Product Registration

Reg. # 18975 issued by Italian Ministry of Health

### Packaging

40 units/box. The box dispenser can be wall-mounted.


6 box/outercase (240 units/outer)

24 outers/pallet (5760 units/pallet)

Dim: 800x1200x2000mm (EPAL)

Weight: 250kg (approx.)

**Product Ref: IODIO P2**



Ergonomic shape to allow a perfect grip  
High flexibility to ensure accurate cleaning  
Soft and kind to the skin  
With nailcleaner for accurate cleaning under the nails  
Advanced bristles system  
Easy peel-off single wrapping  
Single use

**NEX**  
MEDICAL Surgical<sup>®</sup>scrub

[www.nexmedical.it](http://www.nexmedical.it) [info@nexmedical.it](mailto:info@nexmedical.it)

Made in ITALY

# NEX IODIO P2 SURGICAL® SCRUB

Surgical brush/sponge with nailcleaner impregnated with Povidone-Iodine antiseptic solution

## TECHNICAL SPECIFICATION



### 1. PRODUCT DESCRIPTION

**Product Ref:** NEX IODIO P2

**Type of Product:** surgical brush/sponge with nailcleaner impregnated with approx 20ml of Povidone-Iodine antiseptic solution

**Active ingredient(s):** Povidone-Iodine 7,5% (0,75% free Iodine)

**Single use**

**Colour Cod.:** Brown

### 2. PRODUCT REGISTRATION

Reg. # **18975** issued by Italian Ministry of Health  
D.G.F.D.M./405 of 17 Nov. 2004

### 3. INDICATIONS

Surgical scrubbing and washing of hands for pre-operative antiseptic procedures.

### 4. TECHNICAL FEATURES

Dimensions: 45mmx90mmx35mm

Packaging: individually wrapped

### 5. MATERIALS

**Surgical Brush/Sponge:**

Brush: PE (polyethylene) medical grade LDPE (low density polyethylene)

Sponge: Polyurethane

**Nailcleaner:**

PP (polypropylene) medical grade

**Hot melt adhesive:**

medical grade (reg. FDA CFR 21.175.105)

All components are LATEX free.

### 6. FORMULATION COMPOSITION

**Ingredients for 100g:**

Povidone Iodine g, 7,5

Coformulants and water q.s. to g. 100

PH: Min (4,5) ; Max (6,0)

Solids dried @105°C: Min (12,02); Max (14,7)

Sp.Gr.: Min (1,025); Max (1,045)

Titration α: Min (6,75); Max (8,25)

### 7. DERMAL IRRITATION TESTS

Dermal irritation and skin sensitisation:

**OECD guidelines N° 404 and 410 and EN ISO 10993**

**Directive 67/548/EEC ECB. JRC testing method B.6**

### 8. MICROBIOLOGICAL PROPERTIES

**Bactericidal: EN 1040, EN 12054** (Surgical handrub and handwash) and  
**EN 12791** (surgical hand disinfection)

### 9. DIRECTIONS FOR USE

Wet hands and forearms with water. Scrub for at least 3 minutes paying particular attention to the nails, cuticles and interdigital spaces.

### 10. STORAGE INSTRUCTIONS

Avoid storing at excessive heat. Protect from light.

### 11. EXPIRY

3 years

### 12. MANUFACTURER

Nex Medical sas

Settimo Milanese - Milan (ITALY)

registration # PMC N. 20/2008/Off.280 issued by Italian Ministry of Health  
FDA DUNS NO 43-356-4874, LC 50006

### 13. PACKAGING

**Primary pack:** flexible blister pack unit

Materials:

Blister: flexible thermoformed PE medical grade film

Top-foil: flexible PE medical grade film

(FDA 21 CFR: 177/1630-175/105-177/1520)

Features: easy peel-off system

**Secondary pack:** display box-dispenser 40 units

Materials:

Cardboard

Dimensions: 145mm x 125mm x 535mm

Features: front easy opening; inspection window;

eurohole for wall-mounting

**Outercase:** corrugated outercase for 6 box-dispensers

Materials:

Corrugated C flute

Dimensions: 300mm x 400mm x 550mm

Features: american std box

Primary Pack	Secondary Pack	Outercase	Pallet 800x1200mm
Blister single unit	Dispenser (40pcs)	6 dispensers (240 pcs)	24 outers (5760 pcs) h. 200cm

### 14. LABELLING

**Primary pack & box dispenser (printed outercase (labelled):**

product description

product code; disposable logo;

production batch - date - exp;

manufacturer's name and address;

directions for use and washing procedures (on box-dispenser)

**NEX**  
MEDICAL Surgical® scrub

www.nexmedical.it info@nexmedical.it

Made in ITALY

