



Textile gloves

Article No. 281731-281736



Description

Knitted glove in 100% white cotton
Latex free rubber band at the handcuff
White coloured glove with colour coded rubber bands for each size
Unsterile
Avoid storage in direct sunlight.
To be hand washed in a mild soap. Avoid drying in a tumbler.
Produced for Abena A/S, Denmark, in China under ISO 9001

CE

CE-labelled CAT I EU 2016/425 For minimum risk only.

Art. No.	Product	Size	Palm mm	Length mm	Quantity Dispenser/carton
281731	XS - Blue	6	85	200	12/300
281734	S - Red	7	85	210	12/300
281736	M - Orange	8	85	220	12/300
281735	L - White	9	85	240	12/300
281732	XL - Yellow	10	90	245	12/300
281733	XXL - Black	11	90	250	12/300

Storage and handling

Storage: To be stored at a dry and dark place in the original packaging.
Handling: Avoid heat, humidity, bright light and ozone