



# Flocare ENFit adaptere



**293688 / 589732: ENFit → ENLock adapter: til at forbinde ENFit ernæringssæt med ENLock sonde. Inkluderet i Flocare overgangssættet.**

## Fra ENFit ernæringssæt...

293650 Flocare pack set infinity  
293651 Flocare universal set infinity  
293652 Flocare pack set gravity  
293653 Flocare universal set gravity  
293654 Flocare pack mobile set infinity  
293655 Flocare pack set twin line infinity  
293656 Flocare top till reservoir set gravity



## eller fra ENFit sprøjter...

293669 Ernæringssprøjte Flocare 1ml  
293670 Ernæringssprøjte Flocare 3ml  
293671 Ernæringssprøjte Flocare 5ml  
293672 Ernæringssprøjte Flocare 10ml  
293673 Ernæringssprøjte Flocare 20ml  
293674 Ernæringssprøjte Flocare 60ml  
293675 Ernæringssprøjte Flocare 100ml

## ...til ENLock sonder

293550 Flocare pur tube multi ch5/50cm  
293551 Flocare pur tube multi ch6/60cm  
293552 Flocare pur tube multi ch6/110cm  
227357 Flocare pur sonde ch 8/110cm  
227355 Flocare multi purpose ch 10/110cm  
227360 Flocare multi purpose ch 12/110cm  
227356 Flocare multi purpose ch 14/110cm  
293565 Flocare Bengmark ch 08/145cm  
293566 Flocare Bengmark PEG/PEJ ch09  
227353 Flocare Bengmark ch 10/145cm  
293570 Flocare PEG gastrostomisonde ch10  
293571 Flocare PEG gastrostomisonde ch14  
293572 Flocare PEG gastrostomisonde ch18