



Vinyl Gloves, Powder-free

Art. Number 4431, 4424-4426, 4430

Classic

Vinyl gloves, Powder free, Smooth



Description

Examination gloves Vinyl

Powder-free

Non-sterile

Transparent

Smooth surface

With rolled cuff

Plasticizer: DINP

Recommend for handling of non-fatty foodstuff

Produced for Abena A/S, Denmark, in China under ISO 9001 and ISO 13485

CE

Dual CE-labelled according to:

Medical Device regulation EU 2017/745 Class I

PPE Regulation EU2016/425, Type B CAT III

AQL 1.5 according to the requirements in the harmonised standard EN 455-1

Standards

EN 455-1, 2, 3, 4

EN ISO 374-1: 2016 KPT

EN ISO 374-4: 2019

EN ISO 374-5: 2016 VIRUS

EN ISO 21420:2020

EN 16523-1: 2015

Thickness

Palm min. 0,05 mm

Fingers min. 0,07 mm

Physical properties

Force at break

Before challenge 3.6 Newton

After challenge testing at 70°C for 166 hours 3.6 Newton

Elongation

Before challenge min. 300 %

After challenge testing at 70°C for 166 hours min. 300 %



Vinyl Gloves, Powder-free

Art. Number 4431, 4424-4426, 4430

Tests performed

Air Leak Test and Water Leak Test
Sensitization test
Irritation test
Cytotoxicity
Force at break
Powder residuals
Bacteria and Fungi test
Viral penetration
Food approval (please see separate food datasheet)

Composition

Chemical	CAS
PVC	-
DINP	28553-12-0
Viscosity reducer JNJ	64742-47-8
PU-coating	-
Ca-Zn Stabilizer	-

Art. No.	Product	Size	Palm mm (+/- 5)	Length mm	Cartons per pallet	Quantity Dispenser/carton
4431	X-small	5-6	75	240	40	100/1000
4424	Small	6-7	85	240	40	100/1000
4425	Medium	7-8	95	240	40	100/1000
4426	Large	8-9	105	240	40	100/1000
4430	X-Large	9-10	115	240	40	100/1000

Storage and handling

Storage: 10-30°C
Handling: Avoid heat, dampness, bright light and ozone
Shelf life: 5 years from production date