



# ΑΒΧΔΕΦ

OWLERTON GREEN . SHEFFIELD . S6 2BJ  
TELEPHONE: (0114) 234 4231 FAX: (0114) 231 4966  
Email:- QA@SWANN-MORTON.COM

FM 14713

## NON STERILE CARBON SURGICAL BLADE SPECIFICATIONS

### BLADE - CARBON STEEL SURGICAL GRADE

STEEL TYPE	:	CARBON	
SPECIFICATION	:	BS 2982 : 1992	
ANALYSIS	:	CARBON	1.20% - 1.30%
		SILICON	0.10% - 0.35%
		MANGANESE	0.20% - 0.45%
		CHROMIUM	0.10% - 0.40%
		SULPHUR	0.025% max.
		PHOSPHORUS	0.035% max.

### 5 packs - Codes 0101 –0134

FOIL	PRINTED 30 MICRON ALUMINIUM LAMINATE HDPE, 25 MICRONS
RUST INHIBITOR	SODIUM BENSOLATE IMPREGNATED V.P.I. KRAFT PAPER
BOX	600 MICRON COATED WHITE LINED CHIPBOARD

### Single Peel packs - Codes 3001 – 3092

TOP WEB	40gsm PAPER WITH 50 MICRON LOW DENSITY POLYTHENE
BOTTOM WEB	60gsm PAPER WITH 50 MICRON LOW DENSITY POLYTHENE 1g PW213a WITH 6g PW GRID LACQUER
RUST INHIBITOR	SODIUM BENSOLATE IMPREGNATED V.P.I. KRAFT PAPER(BROWN)
<b>OR</b>	
ALTERNATIVE BOTTOM WEB	PET 12 MICRON & ADHESIVE LAQUER 50 MICRON COEXTRUDED POLYPROPYLENE
ALTERNATIVE RUST INHIBITOR	SODIUM BENSOLATE IMPREGNATED BLEACHED V.P.I. (WHITE)